



Certificate of Conformity

VWR LID FOR CRUCIBLES

VWR International certifies that the products comply in every respect with the relevant particular specifications and/or drawings and/or relevant standards and regulations in force, mentioned in this document.

VWR International certifie que les produits sont conformes à tous égards aux spécifications particulières appropriées et / ou dessins, normes et règlements en vigueur, mentionnées dans le présent document

VWR International bestätigt, dass das Produkte allen zutreffenden Spezifikations- und/oder Zeichnungs- und/oder Standardanforderungen sowie den im Dokument aufgeführten gültigen Vorschriften entspricht.

VWR Cat. No.	For Ø (mm)	Pk
459-0215	30	5
459-0216	35	5
459-0217	40	5
459-0218	45	5
459-0219	50	5
459-0220	60	5
459-0221	70	5
459-0222	80	5

Intended uses: Laboratory applications

Material: Porcelain C110

Shape: With handle

Compliant to: ISO material for the manufacture lies within the group C110 according to ČSN EN 60672 –3
Product dimension tolerances are controlled by DIN 40 680 and DIN 12906
Acid resistance acc. to ISO 1775

Properties:

Unglazed products can be used at temperatures of up to 1350°C / 2462°F and glazed ones at temperatures of up to 1000°C / 1,832°F

Products have high resistance to acids and bases

Density 2.4 g/cm³

Certificate of Conformity

Bending strength of unglazed products (70MPa)

Thermal conductivity ($\text{W.m}^{-1}\text{.K}^{-1}$) = 1.3

The compliance of demands is proved by regular tests and followed noticed tests from an independent laboratory:

Test according ČSN EN 60672-2

Test resistance to acids and bases

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **August 2018**