

Certificate of Conformity

VWR Catalogue Number	See List
Description	Sleeve protectors, VWR® Basic
Country of Origin	Manufactured in China
Date of Issue (yyyy-mm-dd)	2022-12-23

20 µm

Catalogue Number	Colour	Length [mm]
113-9073	White	400
113-9074	Blue	400
113-9075	Yellow	400
113-9076	Green	400
113-9077	Red	400

40 µm

Catalogue Number	Colour	Length [mm]
113-9078	Blue	400

Quality System Compliance

Products are manufactured under the **ISO 9001:2015 & ISO 13485:2016** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected according to 4.0 AQL level II sampling plan

For 20 µm

	Test / Standard	MD	CD
Break Tensile Properties [N]	ISO 13934-1	>20	>10
Elongation Rate [%]	ISO 13934-1	>100	>100

For 40 µm

	Test / Standard	MD	CD
Break Tensile Properties [N]	ISO 13934-1	>30	>15
Elongation Rate [%]	ISO 13934-1	>200	>200

Product Specifications

Material	PE Film.
Sterilization	Non-Sterile.
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	Nature rubber latex but without exposure to outside.
BPA Statement	Not Applicable.
DEHP Statement	Not Applicable.
RoHS	Not Applicable.

REACH Statement Not Applicable.

Storage Conditions Temp: 0- 50°C, Humidity: 45-85% without sunshine exposure.

Shelf Life Not Applicable.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoQ using new layout	12-2022