



TECHNICAL DATA SHEET

HYGIENE MASK 3 PLY – Green – earloops 113-0714

CHARACTERISTICS:

Hygiene masks with ties are made of non-woven and melt-blown polypropylene material, commonly used in a clean environment. Face Masks are manufactured and packed in a stringent clean environment.

PRODUCTS DESCRIPTION: Face Mask with tie, glass free, lint-free, low particulating, no exposed metal

PHYSICAL PROPERTIES:

Material:	
Front Layer:	non-woven Polypropylene (green 20 gram/m ²)
Back Layer:	non-woven Polypropylene short fiber (white 25 gram/ m ²)
Center:	Melt-blown Polypropylene (white 20 gram/ m ²)
Ties:	non-woven Polypropylene (white 38 gram/ m ²)
Size:	
Length:	175mm +/- 2mm
Height:	95mm +/- 2mm
Nosepiece:	115mm +/- 2mm
Ties:	420mm +/-5mm
Colors	Green outer, white inner, white ties
Bacterial Filtration Efficiency:	> 90 %
Particle Filtration Efficiency:	> 90 %
Particle size:	0.1Micron

PACKING:

50 pieces/box, 40 boxes/carton

QUALITY ASSURANCE CRITERIA:

Face masks are guaranteed to meet 4.0 AQL level II sampling plan per MIL-STD-105E for the following defects criteria:

- | | | | |
|----|------------------------------------|----|-----------------------|
| 1. | Ultrasonically welded not strong | 2. | Irregular PP material |
| 3. | Dirty PP material | 4. | Nosepiece exposed |
| 5. | Physical dimensions non-compliance | | |

December 6, 2006