

Product Information Package

PRODUCT NAME	D(-)-FRUCTOSE PH.EUR.
PRODUCT CODE	24282
GRADE	PH.EUR
VERSION	0
DATE	04-10-2022

Table of Contents

TABLE OF CONTENTS	2
1. GENERAL INFORMATION	3
2. REGULATORY INFORMATION	3
3. SUPPLY CHAIN	3
4. MANUFACTURER	4
5. FILLING AND PACKAGING	4
6. CERTIFICATE OF ANALYSIS (COA)	6
7. SHELF LIFE	6
8. LABEL	6
9. PRODUCT INFORMATION	7
10. PRODUCT DOCUMENTATION	7
PRODUCT STATEMENTS	8
VWR FILLING AND PACKAGING PROCESS FLOW	11
SUPPLY CHAIN FLOW	12
PACKAGING SPECIFICATIONS	13
PACKAGING SPECIFICATIONS	14
PACKAGING SPECIFICATIONS	15

1. General Information

VWR Part Number	24282
VWR Product Name	D(-)-FRUCTOSE PH.EUR.
CAS number	57-48-7
Origin	Vegetable

2. Regulatory Information

Is the finished product compliant with the following Pharmacopoeia monographs?	
• Ph. Eur.	Yes
• USP	Not tested
• JP	Not tested
• BP	Not tested
• NF	Not tested
• Other (see statement)	N/A
Does VWR hold any regulatory file (CEP, DMF) on the product?	No

3. Supply Chain

With respect to the product in question, does VWR perform:	
• Production	No
• Filling, packaging, relabeling	Yes
• Quality Control	Yes
• Storage and Distribution	Yes
State the name and address of the original manufacturer The original manufacturer of the material can be made available upon request. The name and address will be disclosed subsequent to the signing of a confidentiality agreement with VWR.	
If VWR is not the original manufacturer, is the product delivered to VWR:	
• Directly from manufacturer?	No
• Always from the same manufacturing site?	No
• From an approved supplier	Yes
For a detailed supply chain description see the Supply chain statement in attachment	

4. Manufacturer

Is the manufacturer certified according to	
• ISO 9001	Yes
• ISO 14001	No
• GMP	Yes
• EXCiPACT	No
• Other	N/A
Is the production site dedicated for the product?	
No	
If not, other products?	Citrates, Sweeteners, Gluconates, Phosphates
• Are materials that are high effective substances and/or toxic produced on the same manufacturing site?	No
• Are hormones produced on the same manufacturing site?	No
• Are beta-lactam antibiotics produced on the same manufacturing site?	No
• Are cytostatics produced on the same manufacturing site?	No
Is the production process continuous or batch driven?	
Batch driven	

5. Filling and Packaging

State the VWR Filling and Packaging site:	
VWR International bv - Part of Avantor, Researchpark Haasrode 2020, Geldenaaksebaan 464, B-3001 Leuven, Belgium	
What quality system is applied to this product?	
• ISO 9001	Yes
• ISO 14001	Yes
• ISO 45001	Yes
What level of GMP is applied to this product?	See GMP statement
All certifications of the production site:	
• ISO 9001	Yes
• ISO 14001	Yes

• ISO 45001	Yes
• EU-GMP (limited product scope)	GMP laboratory, limited scope
• ISO17034 (limited product scope)	Yes
Describe the quality management system: The Quality Management System is described in detail in our Self Assesment.	
Is the product filled under controlled conditions?	
• At least Cleanroom ISO Class 8 at rest (0.5µm; 5µm)	Yes
• Laminar flow	No
Is the filling process automated or manual?	Manual
Is the filling room dedicated or multiproduct during filling?	Dedicated
Is a gowning procedure followed to enter the cleanrooms?	Yes
Do you have a system for reconciliation of labels?	Yes
Can a batch be traced back to :	
• The bulk material	Yes
• QC release	Yes
• Filling documents	Yes
How long are batch records retained? Records are retained for at least one year past the expiry date.	
What type of water is used for cleaning purposes? Potable water and final rinse with Purified water compliant to EP and USP	
What is your definition of a batch? A batch is defined as a product volume produced from one manufacturer batch during a fixed period of time	
Lot number system? VWR lot numbers are unique, nine digit alpha-numeric codes compromised of the following: The first two numbes are the last two digits of the year of issuing the filling order. The letter indicates the month (A= Jan, B= Feb, etc.) in which the filling order was issued. The next two numbers indicate the day of printing of the filling order. Finally, the last four digits increase sequentially. This allows to estimate the period in which the product was filled.	

6. Certificate of Analysis (COA)

Which information is stated on the VWR COA	
• VWR Article number	Yes
• Product name	Yes
• Batch number	Yes
• Expiry date	Yes
• Production date	Yes
• Filling date	No
• Specifications	Yes
• Results	Yes
• Electronic signature of QC Manager	Yes

7. Shelf Life

What is the shelf life of the product?	24
How is the expiry date attributed?	copied from supplier
Is the stated shelf life based on	
• Stability testing	Manufacturer information
• Retesting retained samples	No
• Literature	No
• Historical data	No

8. Label

Which information is stated on the VWR product label?	
• VWR product code	Yes
• Product name	Yes
• Batch number	Yes
• Expiry date	Yes
• Production date	No
• Filling date	No
• Address of the VWR filling site	Yes
• Specifications	Yes
• Results	Yes

9. Product Information

We point out that VWR International does not perform these specific testing on the product. But to the best of our knowledge, taking into consideration the raw materials used and/or the manufacturing process, the product is expected to be free from:

Dioxines	No information available
Mycotoxins	No information available
Phtalates	No information available
Melamine	No information available
Aflatoxin	No information available
Nanomaterials	No information available
Asbestos	No information available
Latex	No information available
Carcinogenic, mutagenic and reprotoxic substances	No information available
Pesticide	No information available

10. Product Documentation

Specification	Available on www.VWR.com
SDS	Available on www.VWR.com
BSE/TSE Statement	Available
GMO Statement	Available
Residual solvent Statement	Not available
Residual Metal Catalyst and Elemental Impurities	Available
GMP Statement	Available
Nitrosamine statement	Available
Allergen Statement	Available
VWR Filling and Packaging Process Flow	Available
Manufacturing Process Flow	Not available
Supply chain statement	Available
Packaging specifications	Available
Other	N/A

Product Statements

24282 - D(-)-FRUCTOSE PH.EUR.

BSE/TSE

With reference to guideline EMA/410/01 Rev. 3 of the EC, we declare that this product is manufactured without the use of raw materials of animal or human origin. During processing the product does not come in contact with animal material. Therefore the product is considered safe regarding TSE and BSE.

GMO

With reference to Regulation 1829/2003 and 1830/2003 concerning the labeling and traceability of genetically modified organisms (GMO) and of food and feed produced from genetically modified organisms, we hereby declare that this product is neither produced from, nor consists of genetically modified organisms. Hence, to the best of our knowledge, there is no objective reason to suspect a risk associated with GMO in the product.

Allergen

To the best of our knowledge, taking into consideration the raw materials used and the manufacturing process, allergenic products according to Regulation (EU) No 1169/2011 Annex II are not expected to be present in the final product except for

None ;

For your information, the allergenic products included in Regulation (EU) No 1169/2011 Annex II are:

- Cereals containing gluten and products thereof;
- Crustaceans and products thereof;
- Eggs and products thereof;
- Fish and products thereof;
- Peanuts and products thereof;
- Soybeans and products thereof;
- Milk and products thereof (including lactose);
- Nuts and products thereof;
- Celery and products thereof;

- Mustard and products thereof;
- Sesame seeds and products thereof;
- Sulphur dioxide and sulphites at concentrations of more than 10 mg/kg;
- Lupin and products thereof;

We point out that VWR International does not perform any testing on allergens in the product.

Residual Metal Catalyst and Elemental Impurities

Based on the implementation of ICH Q3D into the European Pharmacopoeia and US Pharmacopoeia, all elements from ICH Q3D Class 1 and 2A (i.e. Cd, Pb, As and Hg plus Co, V and Ni) are listed on our Certificates of Analysis and limits of their maximum ppm concentration will be shown. In addition, the elements of Class 2B or 3 are reported in case they are intentionally added within the manufacturing process. Those elements are analyzed by Atomic Absorption Spectroscopy or Inductively Coupled Plasma using Optical Emission Spectroscopy or Mass Spectroscopy detection.

Based on the results gathered so far, we implemented skip lot testing for the metal catalyst residues. This means we will perform a full analysis on the residues of metal impurities every 5 lots or at least every 2 years if not ordered at high frequency. In case skip lot testing is applied for metal catalyst analysis of a concerned lot, the results of the latest performed analysis will be copied.

This process will not have any impact on the Certificate of analysis and we will continue to report the elements of Class 1 and 2A. Class 2B or 3 are reported in case they are intentionally added within the manufacturing process.

The specifications were set based on our quantitation limits of the ICP method.

The risk of having a non-compliant result reported is mitigated by:

- Qualified source (compendial grade) with ICH Q3D statement
- Enough data have been gathered over the last 5 years to support skip lot testing
- A full analysis is performed every 5 lots or at least every 2 years for less frequent orders.
- In case a metal of Class 1 and 2A is a specification set by the monograph, we will test it and not apply skip lot testing.
- The ICHQ3D guideline does not state that we need to measure the residues of metal catalysts on each lot.

GMP

Our pharma-grade materials are prepared according to principles of good manufacturing practice (GMP). These principles are expected to be applied for the manufacture and repack of products with a designation of a pharmacopoeia in the title of the substance. The applied principles of GMP for pharmaceutical

products are based on IPEC - PQG Joint Good Manufacturing Practices Guide for Pharmaceutical Excipients. The IPEC - PQG GMP Guide is also the basis for our supplier audits. Processes are verified to comply with an adequate quality level of manufacturing and quality control. We trust that the applied GMP principles cover the quality management system and the extent of GMP necessary throughout manufacturing processes for pharma-grade products.

Nitrosamine

Following an internal evaluation in accordance with regulatory authority guidance and as outlined in the ICH Q9 guideline on quality and risk management, VWR International ("VWR") has determined that nitrosamines and related compounds are not likely to be present above reportable limits in this product. VWR's review of the Products for risk of nitrosamine included the assessment of potential nitrogen sources and other chemistries with potential for nitrosamine impurity formation. Following such assessment, in line with the EMA notice "Information on nitrosamines for marketing authorization holders" (EMA/189634/2019) by the European Medicines Agency, dated 19 September 2019, and its related Questions and Answers document (EMA/CHMP/428592/2019 Rev. 2), dated 20 December 2019, VWR is able to confirm as follows:


1. None of our manufacturing processes utilize sodium nitrite or other nitrosating agents as intentionally added ingredients in the presence of secondary, tertiary amines or quaternary ammonium salts within the same or different manufacturing step.
2. None of our manufacturing processes utilize sodium nitrite or other nitrosating agents as intentionally added ingredients in combination with reagents, solvents or catalysts, which are susceptible to degradation to secondary or tertiary amines, within the same or different manufacturing step.
3. During our manufacturing process the products are not coming into contact with nitrosamines or nitrosamines containing products.
4. During our manufacturing process the products are not coming into contact with equipment used to handle nitrosamines or nitrosamines containing products.
5. We do not outsource any nitrogen containing solvents or other nitrogen containing reagents to 3rd parties for recovery or reprocessing.
6. None of our products are sold in blister packaging containing nitrocellulose printer primer that may react with amines in printing ink to potentially generate nitrosamines.

Considering the very low potential for the presence of nitrosamine in VWR manufactured and/or repacked products as set forth above, no specific testing was or is intended to be performed by VWR to confirm the absence of nitrosamine or related compounds in the Products. VWR makes no warranties or representation, express or implied, related to the presence or absence of nitrosamines or related compounds in any of the Products.

VWR Filling and Packaging Process Flow

24282 D(-)-FRUCTOSE PH.EUR.

The flowchart below gives an overview of the filling and packaging flow at VWR International.



Supply chain Flow

24282 D(-)-FRUCTOSE PH.EUR.

With regard to this product, we hereby confirm that the current supply chain of the product is as follows:



Company	Activity
Original Manufacturer	Manufacturing ;QC testing ;Labelling ;Warehousing ;Distribution ;
Supplier	Warehousing ; Distribution ;
VWR International bv Researchpark Haasrode 2020, Geldenaaksebaan 464, B-3001 Leuven Belgium	Repacking ;Relabelling ;QC testing ;
VWR Distribution Warehouses	Warehousing ; Distribution ;

Depending on the customer requirements, various VWR warehouse can be added to the supply chain.

Packaging Specifications

24282 D(-)-FRUCTOSE PH.EUR.

Part number:	24282.290
Packed quantity:	1 KG
Primary Container	
Type	BOTTLE UN X SOLIDS (EC 10/2011 & EC 1935/2004)
Measurements	1L DIN60
Material	HDPE WHITE
Secondary Container	
Type	NA
Material	NA
Primary Closure and Seal	
Type	SCREW CAP UN tamper evident
Measurements	DIN60
Material	PP BLEU
Secondary Closure	
Type	NA
Material	NA

Packaging Specifications

24282 D(-)-FRUCTOSE PH.EUR.

Part number:	24282.368
Packed quantity:	5 KG
Primary Container	
Type	BOTTLE UN X SOLIDS (EC 10/2011 & EC 1935/2004)
Measurements	6L DIN80
Material	HDPE WHITE
Secondary Container	
Type	NA
Material	NA
Primary Closure and Seal	
Type	SCREW CAP WITH.C. (EC 10/2011 & EC 1935/2004)
Measurements	DIN80
Material	PP BLUE
Secondary Closure	
Type	NA
Material	NA

Packaging Specifications

24282 D(-)-FRUCTOSE PH.EUR.

Part number:	24282.461
Packed quantity:	25 KG
Primary Container	
Type	DUSTFREE BAG (EC 10/2011 & EC 1935/2004)
Measurements	900X720MMX100µM
Material	LDPE
Secondary Container	
Type	BUCKET UN Y (EC 10/2011 & EC 1935/2004)
Material	PP WHITE
Primary Closure and Seal	
Type	CABLE TIE (STRAP)
Measurements	3,5X215
Material	PLASTIC
Secondary Closure	
Type	COVER FOR BUCKET UN
Material	BLUE

This document is created based on the information kept in our quality database.

By *QA VWR chemicals - Willem Geerts*

On *October 4, 2022*