

VWR® Cleanroom Paper

01. PROPRIETARY BINDING
TECHNOLOGY

02. SUITABLE FOR USE IN ISO 4 TO 6
(CLASS 10 TO CLASS 1000) CLEANROOM

03. RECYCLABLE



VWR® CLEANROOM PAPER


Cleanroom paper is designed for writing, printing, and photocopying in cleanroom environments where the prevention of chemical exportation and particle contamination is critical.

- Proprietary binding technology
- Suitable for use in ISO 4 to 6 (Class 10 to Class 1000) cleanroom
- Produced in ISO 4 (Class 10) cleanroom
- Excellent printability and heat toner adhesion
- High opacity and durability
- Recyclable
- Autoclavable

Specifications

Weight	80 g/m ² (22 lbs.)	100 g/m ² (28 lbs.)
Basis Weight g/m ²	78-83	102-107
Caliper, Micron	94-104	110-120
Tear Strength, MD/CD	74/75	80/81

Description	Paper Size, in.	Weight, g/m ² (lbs.)	Packaging	Cat. No.	Case of
Synthetic Latex					
White	8.5 x 11	80 (22)	250/10/2500	87004-188	2,500
White	8.5 x 11	105 (28)	250/10/2500	94001-262	2,500
Seafoam	8.5 x 11	105 (28)	250/10/2500	94001-264	2,500
Blue	8.5 x 11	80 (22)	250/10/2500	87004-191	2,500
Green	8.5 x 11	80 (22)	250/10/2500	87004-194	2,500
Pink	8.5 x 11	80 (22)	250/10/2500	89038-762	2,500
Yellow	8.5 x 11	80 (22)	250/10/2500	87004-192	2,500



Description	Paper Size, in.	Weight, g/m ² (lbs.)	Packaging	Cat. No.	Case of
Polyethylene (Non-Latex)					
White	8.5 x 11	80 (22)	250/10/2500	87003-434	2,500
White, 3-Hole	8.5 x 11	80 (22)	250/10/2500	87003-430	2,500
White	8.5 x 11	105 (28)	250/10/2500	87004-170	2,500
White	11 x 17	80 (22)	250/5/1250	89201-140	1,250
White	11 x 17	105 (28)	250/5/1250	89201-142	1,250
Blue	8.5 x 11	80 (22)	250/10/2500	87003-436	2,500
Blue, 3-Hole	8.5 x 11	80 (22)	250/10/2500	87003-432	2,500
Blue	11 x 17	80 (22)	250/5/1250	89201-138	1,250
Cream	8.5 x 11	80 (22)	250/10/2500	82030-866	2,500
Green	8.5 x 11	80 (22)	250/10/2500	87004-164	2,500
Pink	8.5 x 11	80 (22)	250/10/2500	82028-074	2,500
Yellow	8.5 x 11	80 (22)	250/10/2500	82028-072	2,500

VWR® CLEANROOM NOTEBOOKS

Used for writing in cleanroom environments, these spiral-bound notebooks feature latex-free paper in college-rule. Covers are constructed of semitransparent, high-density polyethylene. Notebooks have 50 sheets, 100 pages each, and are recyclable.

- Cleanroom processed and packaged
- Amide and silicone-free
- Proprietary binding technology
- Recyclable
- Autoclavable

Specifications

Weight, Synthetic Latex	80 gsm (22 lbs.)
Basis Weight (g/m ²)	78-83
Caliper, Micron	94-104
Tear Strength, MD/CD	74/75



Ask your VWR sales representative about custom printing!

Color	Paper Size, in.	Weight, g/m ² (lbs.)	Cat. No.	Case of
White	3 x 4	80 (22)	10159-918	40
White	5.5 x 8.5	105 (28)	10159-920	20
White	8.5 x 11	105 (28)	10159-922	10