

Certificate of Conformity

VWR Catalogue Number	129-0728
Description	Plain Blender Bags, 400 mL Capacity, with 300×190 mm Dimension
Lot Number	02042023
Date of Manufacture (yyyy-mm-dd)	2023-02-04
Date of Sterilization (yyyy-mm-dd)	2023-02-10
Expiry Date (yyyy-mm-dd)	2026-02-04
Country of Origin	Manufactured in China
Date of Issue (yyyy-mm-dd)	2026-01-27

Quality System Compliance

This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared their Quality Management System (QMS) complies with [ISO 13485:2016](#).

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Blending Test	:	Bags have to be resistant for 3 minutes blending (Frequency: 4 times/second).
Welding strength	:	> 12N/15 mm.
Tensile Strength	:	> 25MP/10 mm.
Elongation	:	>600%.

Product Specifications

Material 100 % virgin raw material (LDPE + HDPE)

Raw material used to produce these products are in accordance with FDA regulation 21 CFR § 177.1520, for use in contact with food applications

Sterilization	Product labeled as sterile is gamma irradiated and dose released upon ISO11137 recommended practices in effect at the time of validation. Products labeled sterile meet a minimum requirement of 10 ⁻³ SAL with a specified dose range of 10–40 kGy.
ATP Assay	Not applicable.
DNase & RNase Free	Not applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	Not applicable.
Latex Statement	Not applicable.
BPA Statement	Bisphenols are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	Not applicable.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.
REACH Statement	Not applicable.

Storage Conditions Room temperature.

Shelf Life 3 years.

Disclaimer: VWR's Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability, and distribution processes. VWR is not a manufacturer and does not engage in manufacturing, co-manufacturing or processing activities such as design, packaging, labeling, testing or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.

VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use.

This document has been produced electronically and is valid without a signature.