

PDS No. 2251xx	PRODUCT DATA SHEET				Page 1 of 1
Revision 01	PS Multipurpose Container, 150 ml				 greiner bio-one
	Greiner Item-No. 2251xx				
Valid for Item-No.:	225170	225180			

1.	Description / Specification	
1.1	Description	PS Multipurpose Container with metal screw cap 225170: standard, no label 225180: standard, plain label
1.2	Dimensions	Container: External length: 107,0 mm Inner diameter: 47,0 mm Outside diameter: 49,0 mm Container with screw cap: External length: 107,7 mm Outside diameter: 54,7 mm
1.3	Volume	150 ml
1.4	Material / Resin	Container: PS (Polystyrene), free from heavy metal Screw cap: Tin plate, free from heavy metal Label: Paper, free from heavy metal
1.5	Colour	Container: clear Screw cap: silver 225180: Label: white
1.6	Sterilization	No
1.7	Quality Control	- Raw Material-Control: physical testing - Product-Control: testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	- Single use only - Aseptically produced (Class II Clean Room)

2.	Features	
2.1	Basic features	-
2.2	Temperature range	-20°C to +60°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	Container: see homepage: www.gbo.com/bioscience →Products→Literature →Technical Information →Chemical Resistance of Resins
2.6	Shelf life	N/A
2.7	Other Information	-

3.	Packaging	
3.1	Pieces / Bag	20
3.2	Pieces / Box	120
3.3	Lot-No.	YYWWXX (year, week, day, shift)
3.4	Other Information	-

4.	Other Information	
	-	

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision	Date	Date	Date	
-	2 May 2007	4 May 2007	8 May 2007	
Date	Name	Name	Name	
-	S. Kaelberer	Dr. L. Marchetti	A. Schulz	