

TECHNICAL NOTE

VWR[®] Industrial Apparel

VWR[®] INDUSTRIAL APPAREL

- Economical
- First-line barrier
- Maximum breathability

VWR[®] BASIC PROTECTION INDUSTRIAL GARMENTS AND ACCESSORIES

- Economical
- First-line barrier
- Color coding

VWR Basic Protection Industrial garments and accessories are disposable and manufactured from lightweight and breathable fabric. They are an excellent basic barrier with a wide choice of color options to use for color coding everything from mitigating risk in a food safety program to managing critical control points or differentiating various production areas. This product line is ideal for a variety of areas including industrial, food and beverage, production and visitors.






VWR® BASIC PROTECTION INDUSTRIAL COVERALLS

- Tapered collar
- Snap front

Fabric Weight	Size	Color	Packaging	Cat. No.
Industrial Coverall with Zipper Closure and Elastic Cuffs				
20 g/m ²	Small	White	Single Over-bag, 25/CS	75874-030
20 g/m ²	Medium	White	Single Over-bag, 25/CS	75874-032
20 g/m ²	Large	White	Single Over-bag, 25/CS	75874-034
20 g/m ²	X-Large	White	Single Over-bag, 25/CS	75874-036
20 g/m ²	2X-Large	White	Single Over-bag, 25/CS	75874-038
20 g/m ²	3X-Large	White	Single Over-bag, 25/CS	75874-040
20 g/m ²	4X-Large	White	Single Over-bag, 25/CS	75874-042
20 g/m ²	Small	Blue	Single Over-bag, 25/CS	75874-044
20 g/m ²	Medium	Blue	Single Over-bag, 25/CS	75874-254
20 g/m ²	Large	Blue	Single Over-bag, 25/CS	75874-256
20 g/m ²	X-Large	Blue	Single Over-bag, 25/CS	75874-258
20 g/m ²	2X-Large	Blue	Single Over-bag, 25/CS	75874-260
20 g/m ²	3X-Large	Blue	Single Over-bag, 25/CS	75874-262
20 g/m ²	4X-Large	Blue	Single Over-bag, 25/CS	75874-264
20 g/m ²	Small	Green	Single Over-bag, 25/CS	75874-266
20 g/m ²	Medium	Green	Single Over-bag, 25/CS	75874-268
20 g/m ²	Large	Green	Single Over-bag, 25/CS	75874-270
20 g/m ²	X-Large	Green	Single Over-bag, 25/CS	75874-272
20 g/m ²	2X-Large	Green	Single Over-bag, 25/CS	75874-274
20 g/m ²	3X-Large	Green	Single Over-bag, 25/CS	75874-276
20 g/m ²	4X-Large	Green	Single Over-bag, 25/CS	75874-278
20 g/m ²	Small	Red	Single Over-bag, 25/CS	75874-280
20 g/m ²	Medium	Red	Single Over-bag, 25/CS	75874-282
20 g/m ²	Large	Red	Single Over-bag, 25/CS	75874-284
20 g/m ²	X-Large	Red	Single Over-bag, 25/CS	75874-286
20 g/m ²	2X-Large	Red	Single Over-bag, 25/CS	75874-288
20 g/m ²	3X-Large	Red	Single Over-bag, 25/CS	75874-290
20 g/m ²	4X-Large	Red	Single Over-bag, 25/CS	75874-292

COVERALL DIMENSIONS

Size	M	L	XL	2X	3X	4X	5X
Body (1)	37½"	38¾"	40"	40½"	43"	43"	43½"
Chest (2)	24½"	25"	27"	29¾"	32"	32½"	33½"
Sleeve (3)	34"	34"	36½"	37"	38"	41"	42"
Leg (4)	28½"	29½"	30"	31"	32"	32½"	33"



VWR® BASIC PROTECTION INDUSTRIAL LAB COATS

- Tapered collar
- Snap front

Fabric Weight	Size	Colour	Packaging	Cat. No.
Lab Coat with Open Cuff and No Pockets				
20 g/m ²	Small	White	Single Over-Bag, 30/CS	76008-420
20 g/m ²	Medium	White	Single Over-Bag, 30/CS	76008-422
20 g/m ²	Large	White	Single Over-Bag, 30/CS	76008-426
20 g/m ²	X-Large	White	Single Over-Bag, 30/CS	76008-428
20 g/m ²	2X-Large	White	Single Over-Bag, 30/CS	76008-430
20 g/m ²	3X-Large	White	Single Over-Bag, 30/CS	76008-432
20 g/m ²	4X-Large	White	Single Over-Bag, 30/CS	76008-434
20 g/m ²	Small	Blue	Single Over-Bag, 30/CS	76008-436
20 g/m ²	Medium	Blue	Single Over-Bag, 30/CS	76008-438
20 g/m ²	Large	Blue	Single Over-Bag, 30/CS	76008-440
20 g/m ²	X-Large	Blue	Single Over-Bag, 30/CS	76008-442
20 g/m ²	2X-Large	Blue	Single Over-Bag, 30/CS	76008-480
20 g/m ²	3X-Large	Blue	Single Over-Bag, 30/CS	76008-482
20 g/m ²	4X-Large	Blue	Single Over-Bag, 30/CS	76008-484
20 g/m ²	Small	Green	Single Over-Bag, 30/CS	76008-486
20 g/m ²	Medium	Green	Single Over-Bag, 30/CS	76008-488
20 g/m ²	Large	Green	Single Over-Bag, 30/CS	76008-490
20 g/m ²	X-Large	Green	Single Over-Bag, 30/CS	76008-492
20 g/m ²	2X-Large	Green	Single Over-Bag, 30/CS	76008-494
20 g/m ²	3X-Large	Green	Single Over-Bag, 30/CS	76008-496
20 g/m ²	4X-Large	Green	Single Over-Bag, 30/CS	76008-498
20 g/m ²	Small	Red	Single Over-Bag, 30/CS	76008-500
20 g/m ²	Medium	Red	Single Over-Bag, 30/CS	76008-502
20 g/m ²	Large	Red	Single Over-Bag, 30/CS	76008-504
20 g/m ²	X-Large	Red	Single Over-Bag, 30/CS	76008-506
20 g/m ²	2X-Large	Red	Single Over-Bag, 30/CS	76008-508
20 g/m ²	3X-Large	Red	Single Over-Bag, 30/CS	76008-510
20 g/m ²	4X-Large	Red	Single Over-Bag, 30/CS	76008-512

MATERIAL SPECIFICATIONS: 20 G SPP

Item	Test method	Standard
Weight (g/m²)	GB/T4669-1995	20.0 ± 1.0
Break tensile properties (N)	MD	ASTM D5034 >25
	CD	ASTM D5034 >16
Elongation rate (%)	MD	ASTM D5034 >30
	CD	ASTM D5034 >30

Fabric Weight	Size	Colour	Packaging	Cat. No.
Lab Coat with Elastic Cuff and Three Pockets				
20 g/m ²	Small	Blue	Single Over-Bag, 30/CS	76218-722
20 g/m ²	Medium	Blue	Single Over-Bag, 30/CS	76218-724
20 g/m ²	Large	Blue	Single Over-Bag, 30/CS	76218-726
20 g/m ²	X-Large	Blue	Single Over-Bag, 30/CS	76218-728
20 g/m ²	2X-Large	Blue	Single Over-Bag, 30/CS	76218-716
20 g/m ²	3X-Large	Blue	Single Over-Bag, 30/CS	76218-718
20 g/m ²	4X-Large	Blue	Single Over-Bag, 30/CS	76218-720
20 g/m ²	Small	White	Single Over-Bag, 30/CS	76204-788
20 g/m ²	Medium	White	Single Over-Bag, 30/CS	76204-790
20 g/m ²	Large	White	Single Over-Bag, 30/CS	76204-792
20 g/m ²	X-Large	White	Single Over-Bag, 30/CS	76204-632
20 g/m ²	2X-Large	White	Single Over-Bag, 30/CS	76204-634
20 g/m ²	3X-Large	White	Single Over-Bag, 30/CS	76204-636
20 g/m ²	4X-Large	White	Single Over-Bag, 30/CS	76204-638
Lab Coat with Elastic Cuff and No Pockets				
20 g/m ²	Small	Blue	Single Over-Bag, 30/CS	76360-482
20 g/m ²	Medium	Blue	Single Over-Bag, 30/CS	76360-484
20 g/m ²	Large	Blue	Single Over-Bag, 30/CS	76360-486
20 g/m ²	X-Large	Blue	Single Over-Bag, 30/CS	76360-488
20 g/m ²	2X-Large	Blue	Single Over-Bag, 30/CS	76360-490
20 g/m ²	3X-Large	Blue	Single Over-Bag, 30/CS	76360-492
20 g/m ²	Small	White	Single Over-Bag, 30/CS	76360-470
20 g/m ²	Medium	White	Single Over-Bag, 30/CS	76360-472
20 g/m ²	Large	White	Single Over-Bag, 30/CS	76360-474
20 g/m ²	X-Large	White	Single Over-Bag, 30/CS	76360-476
20 g/m ²	2X-Large	White	Single Over-Bag, 30/CS	76360-478
20 g/m ²	3X-Large	White	Single Over-Bag, 30/CS	76360-480

VWR® BASIC PROTECTION INDUSTRIAL LAB COATS

Size	S	M	L	XL	2X	3X	4X
Body (in)	33½"	35½"	37½"	40"	43½"	45½"	47½"
Chest (in)	24½"	25½"	26"	26¾"	27½"	28¾"	29"
Sleeve (in)	22"	22"	22"	23"	23"	23"	23"



VWR® BASIC INDUSTRIAL SPP BOUFFANT CAPS

Fabric Weight	Size	Color	Packaging	Cat. No.
Latex-Free				
10 g/m ²	Small	White	100/BG, 10 BG/CS, 1,000/CS	76008-514
10 g/m ²	Large	White	100/BG, 10 BG/CS, 1,000/CS	76008-516
10 g/m ²	Small	Blue	100/BG, 10 BG/CS, 1,000/CS	76008-518
10 g/m ²	Large	Blue	100/BG, 10 BG/CS, 1,000/CS	76008-520
10 g/m ²	Small	Green	100/BG, 10 BG/CS, 1,000/CS	76008-530
10 g/m ²	Large	Green	100/BG, 10 BG/CS, 1,000/CS	76008-532
10 g/m ²	Small	Red	100/BG, 10 BG/CS, 1,000/CS	76008-522
10 g/m ²	Large	Red	100/BG, 10 BG/CS, 1,000/CS	76008-524
10 g/m ²	Small	Yellow	100/BG, 10 BG/CS, 1,000/CS	76008-526
10 g/m ²	Large	Yellow	100/BG, 10 BG/CS, 1,000/CS	76008-528
Pleated Latex-Free				
10 g/m ²	Small	White	100/BG, 10 BG/CS, 1,000/CS	76008-534
10 g/m ²	Large	White	100/BG, 10 BG/CS, 1,000/CS	76008-536
10 g/m ²	Small	Blue	100/BG, 10 BG/CS, 1,000/CS	76008-538
10 g/m ²	Large	Blue	100/BG, 10 BG/CS, 1,000/CS	76008-540
10 g/m ²	Small	Green	100/BG, 10 BG/CS, 1,000/CS	76008-550
10 g/m ²	Large	Green	100/BG, 10 BG/CS, 1,000/CS	76008-552
10 g/m ²	Small	Red	100/BG, 10 BG/CS, 1,000/CS	76008-542
10 g/m ²	Large	Red	100/BG, 10 BG/CS, 1,000/CS	76008-544
10 g/m ²	Small	Yellow	100/BG, 10 BG/CS, 1,000/CS	76008-546
10 g/m ²	Large	Yellow	100/BG, 10 BG/CS, 1,000/CS	76008-548



VWR® BASIC INDUSTRIAL NYLON BOUFFANT CAPS

Fabric Weight	Size	Color	Packaging	Cat. No.
10 g/m ²	X-Small	White	Single Over-Bag, 1,000/CS	75997-646
10 g/m ²	Small	White	Single Over-Bag, 1,000/CS	75829-204
10 g/m ²	Large	White	Single Over-Bag, 1,000/CS	75997-636
10 g/m ²	X-Small	Black	Single Over-Bag, 1,000/CS	76076-784
10 g/m ²	Small	Black	Single Over-Bag, 1,000/CS	76076-782
10 g/m ²	Large	Black	Single Over-Bag, 1,000/CS	76076-786
10 g/m ²	X-Small	Blue	Single Over-Bag, 1,000/CS	76076-778
10 g/m ²	Small	Blue	Single Over-Bag, 1,000/CS	76076-776
10 g/m ²	Large	Blue	Single Over-Bag, 1,000/CS	76076-780

Bouffant Cap dimensions

Size	S	L	XL
Diameter, Relaxed (1)	8"	8½"	9"
Diameter, Stretched (1)	21"	24"	28"



VWR® BASIC PROTECTION INDUSTRIAL BEARD COVERS

Fabric Weight	Size	Color	Packaging	Catalog Number
SPP				
11 g/m ²	Universal	White	100/BG, 10 BG/CS, 1,000/CS	89107-778
11 g/m ²	Universal	White	100/BG, 10 BG/CS, 1,000/CS	89218-950
Nylon				
10 g/m ²	Universal	White	100/BG, 10 BG/CS, 1,000/CS	75829-206

Beard Cover Dimensions

Size	Universal	XL
Diameter, Relaxed (1)	8½"	9"
Diameter, Stretched (1)	18"	18"

ALL TEST RESULTS PROVIDED BY INDEPENDENT THIRD-PARTY TESTING LABORATORIES LOCATED IN THE U.S.

WARNING:

These garments and associated materials are not suitable for use in some environments containing chemicals and/or hazardous agents. It is the responsibility of the user to determine the level of risk in a particular environment and the proper personal protection equipment needed. Garments manufactured from synthetic non-woven material may generate static electricity. Garments that contain an anti-static treatment are not intended to be used as a safety feature. These garments are not recommended to be used in a flammable and/or explosive environment. Contact Avantor® for garment/fabric safety data. The application of these products is out of Avantor's control. Therefore, Avantor makes no warranties, expressed or implied, and assumes no liability as to the performance of these products for a particular use.

CAUTION:

Avoid heat and/or open flame.

TO ORDER

Visit www.avantorsciences.com or call your Avantor Sales Representative today.