

Certificate of Conformity

VWR Catalogue Number	See List
Description	Adhesive PCR film seals, Transparent
Country of Origin	Manufactured in USA
Date of Issue (yyyy-mm-dd)	2022-04-20

Catalogue Number	Description	Sterile	Thickness [µm]
391-1254	Adhesive PCR film seals	Non-Sterile	74.9
391-1255	Adhesive PCR film seals	Sterile	74.9
731-0321	Adhesive PCR film seals to seal 2x8 wells	Non-Sterile	74.9
391-1294	Adhesive PCR film seals	Non-Sterile	92.7

Quality System Compliance

Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material	Polypropylene.
Recommended Temperature	-40 °C to 120 °C.
Sterilization	Product labelled as sterile is gamma irradiated and meets SAL 10 ⁻³ with a required dose range of 18.5-12 kGy. (which is not applicable for Non-sterile products)
ATP Assay	Not applicable.
DNase & RNase Free	This product is free of any detectable DNase/RNase contamination.

DNA Contamination	This product is free of detectable levels of Human/Mouse DNA contamination (PCR Certified).
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	Not applicable.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Not applicable.
DEHP Statement	Not applicable.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr ⁶⁺), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list published 2022-01-17 are used in manufacturing of the raw materials used to produce this product. No routinely analyse is performed.
Storage Conditions	21.1 °C ± 10°; 50% RH ± 10%
Shelf Life	3 years.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.