

Product Specification

Test	Specification
Assay (by GC)	$\geq 95 \%$
Residue after Evaporation	$\leq 2 \text{ ppm}$
Substances Darkened by H ₂ SO ₄	passes test
Water (H ₂ O)	$\leq 0.01 \%$
Ultraviolet Absorbance,(1.00-cm path vs water; curve smooth throughout stated range with no extraneous impurity peaks).	
at 210 nm	≤ 0.2
at 220 nm	≤ 0.08
at 254 nm	≤ 0.008
at 280 nm	≤ 0.005
at 350 nm	≤ 0.005
UV Cut-off, nm	≤ 192

For Laboratory, Research or Manufacturing use.
For use in High Performance Liquid Chromatography.
Filtered through a 0.2 micron filter.Packaged under Nitrogen.



Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000