

Certificate of Conformity

VWR Catalogue Number	See list
Description	VWR®, Micro Cuvettes, Quartz
Country of Origin	Manufactured in China
Date of Issue (yyyy-mm-dd)	2025-11-10

Catalogue Number	Description	Pk
634-0688	Cell, quartz, micro cell, 100 µl, 10 mm path length, beam height 15 mm, 1 pcs	1
634-0689	Cell, quartz, micro cell, 200 µl, 10 mm path length, beam height 15 mm, 1 pcs	1
634-6025	Cell, quartz, micro cell, 500 µl, 10 mm path length, 1 pcs	1

Quality System Compliance

This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared their Quality Management System (QMS) complies with **ISO 9001**.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material	ES Quartz Glass.
Sterilization	Not Applicable.
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.
BSE/TSE	Not Applicable.

Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	Not Applicable.
BPA Statement	Not Applicable.
DEHP Statement	Not Applicable.
RoHS	Not Applicable.
REACH Statement	Not Applicable.
Storage Conditions	The ideal storage environment for cuvettes is a dry, clean, light-protected, and room-temperature area. They should be kept away from moisture (to prevent contamination or mildew), strong light (to avoid affecting their optical performance), and corrosive substances.
Shelf Life	Not Applicable

Disclaimer: VWR's Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability, and distribution processes. VWR is not a manufacturer and does not engage in manufacturing, co-manufacturing or processing activities such as design, packaging, labeling, testing or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.

VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use.

This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoQ using new layout	11-2025