

Product Specification Sheet

VWR IMMERSION THERMOSTAT, MX

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

Cat. No.	Temp. range	Temp. stability	Heating capacity	WxDxH (mm)	Weight	Pk
462-0205	RT*...135 °C	±0.07 °C	1.1 kW	109×97×358	4.5 kg	1

Other characteristics:

- The MX immersion thermostat is user friendly with three control buttons and on-screen prompts.
- The convenient slide control enables the flow rate to be easily adjusted, and it can clamp securely to straight or curved tank walls.
- Large EasyView™ LCD display with alarm or fault icons and English prompts and menus
- User adjustable high and low temperature limits and low liquid level alarms
- Single speed pressure pump
- Integral pump and heater coil protection
- Single point calibration
- Complies with DIN 12876-1, Safety class I

Disclaimer: VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there is any doubt ask us.

Date of issue: **June 2013**