



ZPACK13951065 PLASTIC JERRYCAN LIQ 5 L DIN 42 UN



SPECIFICATIONS

UN-certification :	3H1/Y1,4/100
Material :	HDPE
Neck Type :	DIN 42
Filling Opening :	42,00 mm
Height :	285,00 mm
Length :	187,00 mm
Wide :	137,00 mm

AVAILABLE ACCESSORIES

Adaptator for dispenser DIN 42 :	BRND704349
Classical tap :	29048.098
Fluorinated Tap :	29048.099

Remark : Packaging materials and specifications shown may be modified without preceding notice
Version dd. 22 mai 2012

tel/tél +32(0)16 385 011
fax +32(0)16 385 385
info@be.vwr.com
www.vwr.com

Researchpark Haasrode 2020
Geldenaaksebaan 464 B-3001 Leuven
RPR/RPM Leuven
BTW/TVA BE 0403.593.343