

Calbiochem®

Certificate of Analysis

Protein A, *Staphylococcus aureus*

Batch Number:	D00118510
Material Number:	539202-5MG

Analytical Data

Test	Tolerance	Result
Solubility:		H ₂ O
Appearance:		Lyophilized
Form:		Solid
Purity:	90.0 % - 105.0 %	94.0 %
Binding capacity:	95.0 % - 100.0 %	99.0 %
A320 nm:	≤0.04	0.01

Storage and Handling:	-20°C
-----------------------	-------



Eduardo A. Uribe, Quality Assurance Manager

18-Aug-2011

Date

USA and Canada
Tel (800) 628-8470
technical@calbiochem.com

France
Freephone
0800 126 461

Germany
Freecall
0800 100 3496

Ireland
Toll Free
1800 409 445

Europe
United Kingdom
Freephone
0800 622 935

All Other
European Countries
+44 115 943 0840

All Other Countries
Contact Your Local
Distributor
www.calbiochem.com
technical@calbiochem.com

techservice@merckbio.eu

A Brand of EMD Chemicals Inc., an affiliate of Merck KGaA, Darmstadt, Germany
www.calbiochem.com

FOR RESEARCH USE ONLY. NOT FOR HUMAN OR DIAGNOSTIC USE.