

VWR® BASIC CPE Shoe Covers


VWR Basic CPE Shoe Covers

- Ideal for Short-Wear Applications
- All Film with Heat-Sealed Seams
- Excellent Fluid Resistance

VWR Basic Protection CPE Shoe Covers are manufactured from 8mil thick, low-density chlorinated polyethylene (CPE) fabric, which protects against dirt, grime, and liquid splash. These nonconductive shoe covers have a seamless bottom to repel water. Heat-sealed seams provide excellent fluid resistance. Shoe covers are rigorously tested and manufactured in an ISO Certified facility under stringent process controls to ensure that each product meets exacting quality standards and performs to specification. Our products are validated through independent lab testing.

CHARACTERISTICS:

- Made of 48g CPE film
- Extra large size
- Colour blue or white



VWR Cat.number	Size	A	B	C	Color
113-1491	Extra large	40,64	14,60	8,25	White
113-1492	Extra large	40,64	14,60	8,25	blue

PRODUCT SPECIFICATIONS

MATERIAL PROPERTIES

		Result	Test Standard
PHYSICAL PROPERTIES			
	Weight (g/m2)	48.82	ASTM D3776
	Mullen Bursting Strength (PSI)	6.2	ASTM D3786
	Trapezoid Tearing Strength (Avg. lbs./in.)		
	Lengthwise Yarns	2.4	ASTM D4533
	Widthwise Yarns	3.1	

QUALITY ASSURANCE CRITERIA:

AQL 1,5 per General Inspection Plan level II, single sampling (ANSI/ASQC Z1.4-1993)

- Extraneous material
- Dirty products
- Thread ends loose or too long
- Holes or tears
- Foul odor
- Off colour
- Uneven material thickness
- Hard lump
- Fibre spinning-hook off
- Wrong collar tag
- Physical dimension Non-compliance
- Stick is not correct
- Cloth surface applied wrongly

QUALITY ASSURANCE CRITERIA FOR MAJOR DEFECTS:

AQL 2.5 per General Inspection Plan Level II, single sampling (ANSI/ASQC Z1.4-1993)

- Sewing sundries
- Elastic band end exposed
- Sewing end falling-off
- Jumping stitch
- Missing stitch
- Less than 8 stitches per inch
- Unwanted folds
- Elastic tightness non-consistent