

Product Specification

Material Number: 85704.601
Material description: Phosphate test strips, Dosatest
Date specification set: January 21, 2025

(All batches of this product after this date are tested to this specification. In the unlikely event that the specification is changed a new certificate will be issued)

Gradations

0 - 3 - 10 - 25 - 50 - 100mg/L Phosphate

Color change

Phosphate - white → blue- green

Brand

QUANTOFIX

Additional reagent

Yes

Platform

Semi-quantitative test strips

CE certified

No

Method

Phosphate - Phosphomolybdenum blue

Shelf Life

2.5 years.

Storage Temperature

4–30 °C / 39–86 °F.

Hazardous material

No Information available

Gross weight (incl. packaging)

230.7 g / 0.51 lbs

Product Specification

Packaging dimensions

170 x 124 x 44 mm / 6.69 x 4.88 x 1.73 Inch

Signature

This document has been produced electronically and is valid without a signature.

VWR International
