



## DATA SHEET



**Article-No.** 80.066.906  
**Description:** Dissecting Set

6 parts filled, in a map

**Content:**

- 1 x 80.037.150 - powder spatula, length 150 mm, spatula 9 mm wide, material 1.4301
- 1 x 80.101.100 – microscopy scissors straight, length 100 mm, material carbon steel C50 nickel plated
- 1 x 80.021.150 – spoon spatula, length 150 mm, material 1.4301
- 1 x 80.066.620 – scalpel pointed, length 125 mm, material 1.4021
- 1 x 80.083.010 – section lifter flexible, length 150 mm, material handle 1.4301, blade 1.4310
- 1 x 80.116.130 – tweezers pointed, length 130 mm, material 1.4301