

# Certificate of Conformity

|                                         |                                                                                              |
|-----------------------------------------|----------------------------------------------------------------------------------------------|
| <b>VWR Catalogue Number</b>             | 612-6701                                                                                     |
| <b>Description</b>                      | VWR®, Measuring Flasks with Single Certificate, 5 mL Capacity, 0.040± mL Tolerance, 10/19 NS |
| <b>Lot Number</b>                       | GL2523                                                                                       |
| <b>Date of Manufacture (yyyy-mm-dd)</b> | 2023-06-12                                                                                   |
| <b>Country of Origin</b>                | Manufactured in India                                                                        |
| <b>Date of Issue (yyyy-mm-dd)</b>       | 2025-06-17                                                                                   |

## Quality System Compliance

Products are manufactured in the **ISO 9001:2015** certified plant & calibrated in **ISO 17025:2017** certified laboratory. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

## QC Testing

The product is single inspected and single calibrated in accordance with current applicable QC requirements.

## Product Specifications

The product is in conformance with the requirements of the following standards and regulations.  
Technical standards for products: EN ISO 1042  
Calibration procedure: ISO 4787

|                               |                                                                                                                                 |
|-------------------------------|---------------------------------------------------------------------------------------------------------------------------------|
| <b>Material</b>               | The products are made from DIN ISO 3585:1998, TYPE-1, CLASS- A, Boro 3.3 Glass that meet ASTM E-438, TYPE-1, CLASS- A Standard. |
| <b>Sterilization</b>          | Hot air sterilization, in the oven up to 140 °C<br>Steam sterilization, in an autoclave 121 °C/ 20 min<br>134 °C/ 10 min.       |
| <b>ATP Assay</b>              | Not Applicable.                                                                                                                 |
| <b>DNase &amp; RNase Free</b> | Not Applicable.                                                                                                                 |

|                                |                                                                                                                                                                                                                                                                                                                                                             |
|--------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>BSE/TSE</b>                 | No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE). |
| <b>Cytotoxicity</b>            | Not Applicable.                                                                                                                                                                                                                                                                                                                                             |
| <b>Non-Pyrogenic Statement</b> | Not Applicable.                                                                                                                                                                                                                                                                                                                                             |
| <b>Latex Statement</b>         | Not Applicable.                                                                                                                                                                                                                                                                                                                                             |
| <b>BPA Statement</b>           | Not Applicable.                                                                                                                                                                                                                                                                                                                                             |
| <b>DEHP Statement</b>          | Not Applicable.                                                                                                                                                                                                                                                                                                                                             |
| <b>RoHS</b>                    | Not Applicable.                                                                                                                                                                                                                                                                                                                                             |
| <b>REACH Statement</b>         | Not Applicable.                                                                                                                                                                                                                                                                                                                                             |
| <b>Storage Conditions</b>      | Room temperature.                                                                                                                                                                                                                                                                                                                                           |
| <b>Shelf Life</b>              | No predefined shelf life.                                                                                                                                                                                                                                                                                                                                   |

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.  
This document has been produced electronically and is valid without a signature.