



Sol - M[®] Luer Lock Syringe Trays

- Syringes are packaged with tips facing in one direction for enhanced aseptic technique and workflow automation
- Help improve workflow by reducing the number of single-syringe packaging that needs to be opened
- Available with Sol-M Luer Lock syringes sized 1 through 60mL
- Clear, bold graduations for accurate medication preparation and administration

***samples available upon request**

VWR Cat No.	Product Description	PCS/Case	PCS/Tray
76582-628	Sol-M 1mL Luer Lock Syringe Sterile Convenience Tray	1000	25
76582-798	Sol-M 3mL Luer Lock Syringe Sterile Convenience Tray	300	25
76582-804	Sol-M 5mL Luer Lock Syringe Sterile Convenience Tray	300	25
76582-810	Sol-M 10mL Luer Lock Syringe Sterile Convenience Tray	240	20
76582-823	Sol-M 20mL Luer Lock Syringe Sterile Convenience Tray	120	10
76582-826	Sol-M 30mL Luer Lock Syringe Sterile Convenience Tray	120	10
76582-830	Sol-M 60mL Luer Lock Syringe Sterile Convenience Tray	120	20