


CH0993NG



Eisco Labs NEXTGEN Bunsen Burner, Anti-Tip Design with Handle, with Flame Stabilizer and Gas Adjustment, Natural Gas

(A) Input outside Diameter	(B) Base size	(C) Total Height	(D) Adjustment Type	(E) Additional Feature
0.315" (8 mm)	4.25"x 4.25" (108 x 108 mm)	6.4375" (163.51 mm)	Air flow adjustment sleeve / Brass gas adjustment	Candle stick easy carry handle, Brass Flame Stabilizer