


Dino-Lite software - DinoCapture and DinoXcope

A professional, reliable software environment is essential when working with computer equipment like an USB microscope. All Dino-Lite USB products are delivered with an in-house developed software program.

The Dino-Lite software is continuously developed, is free of charge for Dino-Lite users and has an automatic update feature. Dino-Lite software is available for Windows and for MacOS (DinoXcope). DinoCapture software is intuitive, user-friendly and can be used with hardly any training. Free online & e-mail software support is available.

The DinoCapture 2.0 software is available in many languages, such as: English, German, French, Spanish, Chinese, Japanese, Korean, Portuguese, Italian, Russian, Dutch, Greek, Hungarian, Polish, Romanian, Swedish, Finnish, Danish, Czech, Croatian, Norwegian, Turkish, Arabic.



main software features:

- ▶ capturing photographs, videos or time-lapsed videos
- ▶ saving pictures in several formats
- ▶ advanced image processing
- ▶ measurement options like: line, radius, circle, 3-point circle, angle, etc.
- ▶ measurements on captured images or on live images
- ▶ calibration
- ▶ e-mail integration
- ▶ adding notes and markings on images
- ▶ skype integration for real-time on-line sharing with suppliers, customers or colleagues
- ▶ connect multiple Dino-Lite microscopes
- ▶ controlling lighting options from the software
- ▶ IP functionality for remote viewing of microscopic images
- ▶ barcode/QR code recognition functionality
- ▶ GPS integration

Dino-Lite software

Extended Depth of Field (EDOF)

Extended Dynamic Range (EDR)

Advanced measurement options like:
line, radius, circle, 3-point circle, angle, etc.

Calibration functionality

Flexible LED Control (FLC)

Illumination settings

