



Product Specification Sheet

VWR INOCULATING LOOPS, RIGID, 10 μ l, 500pk

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

VWR Catalog No. : **612-9353**

Length of the Loop: **200 mm $\pm$ 2 mm**
Stick Diameter: **Handle Part: 4 mm $\pm$ 0.6 mm**
Loop Diameter: **Max 3.90 mm**

Top Sheet Material : **Medical Grade Paper Plus 60 GR**
Bottom Sheet Material : **PA/PE Film - 100 μ**

Packaging: **Paper-Film Pouch**
Inner Box Qty: **500**
Outer Box Qty: **2000**

Sterilization Method : **Gamma Irradiation**

Storage Temperature : **2/30°C – Do Not Freeze**

Shelf life : **36 Months**

Disclaimer: *VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there is any doubt ask us.*