

CERTIFICATE OF ANALYSIS



Lot # 0G049

Date of Manufacture July 2020

Recommended Retest date is two years from the date of manufacture

SIZE	1L (34oz)	5L (1.32gal)	20L (5.28gal)	210L (55gal)
DECON #	1002	1003	1004	1005

TEST	SPECIFICATION	UNIT	TEST RESULT
pH, neat solution	12.7 - 13.8	N/A	13.3
pH, 5% v/v solution	11.5 - 12.7	N/A	12.5
Appearance	White turbid liquid	N/A	Pass
Density	Report Result	g/ml	1.090
mEq acid/ gram Contrad 70	0.36 – 0.55	mEq/g	0.49

Reviewed and Approved by:

Date: July 2020



Karen Wend
 Manager, Quality and Technical Services

BSE/TSE Statement: There are no materials of animal origin in the raw materials or in the process aides used to manufacture any Decon products.