

Certificate

Extraction thimble MN 645

REF: 645026

LOT: CF253002

Product description

standard grade made from pure cellulose

typical data	unit	specifies
inner diameter	mm	33
height	mm	205

Extraction thimbles with inner diameter x height of 22x80, 33x94, 33x205, 48x230 mm are in accordance to DIN 12449 for extractors with defined volumes according to DIN 12602 and DIN 12604.

Confirmation

Hereby we confirm that the above-mentioned product has successfully passed our quality control system in accordance with EN ISO 9001 and meets the specific quality criteria.

This document has been produced electronically and is valid without a signature.