



ACROS ORGANICS part of Thermo Fisher Scientific

ACRÖS
ORGANICS

Product Specifications

Product 128191000
1-Naphthol,99+%

General Product Data

Version	2
CAS No	90-15-3
Molecular weight	144.17
Molecular formula	C ₁₀ H ₈ O
Linear formula	C ₁₀ H ₇ OH
Flash point (°C)	125

Product Specifications

Appearance (Color)	Beige to grey to brown
Appearance (Form)	Crystalline powder or crystals or flakes
Infrared spectrum	Conforms
Melting point	95°C to 97°C
GC	>=99.0 %

Product Specifications are subject to amendment and may change over time.

Acros Organics

ENA23, zone 1, nr 1350, Janssen Pharmaceuticaaan 3a, B-2440 Geel, Belgium
Tel +32 14/57.52.11 - Fax +32 14/59.34.34 Internet: <http://www.acros.com>
1 Reagent Lane, Fair Lawn, NJ 07410,USA
Fax 201-796-1329