

PDS No. 103101	PRODUCT DATA SHEET	Page 1 of 1
Revision 04	Tube, 2 ml, PS, Round Bottom	 greiner bio-one
	Item-No. 103101	

1.	Description / Specification	
1.1	Description	PS tube with round bottom, 2 ml Appropriate grip stopper, Cat.-No. 301321
1.2	Dimensions	See Customer drawing
1.3	Volume	<u>Working volume</u> : 1.5 ml
1.4	Material / Resin	PS (Polystyrene), free of heavy metal
1.5	Colour	Clear
1.6	Sterilization	No
1.7	Quality Control	- <u>Raw Material-Control</u> : physical testing - <u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	For single use only

2.	Features	
2.1	Basic features	-
2.2	Temperature range	-20°C to +60°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	5800 x g: swinging-bucket rotor
2.5	Chemical Resistance	See homepage: www.gbo.com/bioscience →Products →Literature →Technical Information→Chemical Resistance of Resins
2.6	Shelf life	N/A
2.7	Other Information	-

3.	Packaging	
3.1	Pieces / Bag	500
3.2	Pieces / Box	3000
3.3	Lot-No.	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	-

4.	Other Information	
	-	

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 03	Date 31 July 2014	Date 1 August 2014	Date 1 August 2014	
Date 17.12.2009	Name S. Kaelberer	Name A. Fuhrmann	Name A. Schulz	