



A

USP VI


8060 NALGENE 50 Platinum-Cured Silicone Tubing

This flexible (50, Shore A), durable, high-purity tubing is designed for a variety of pump and transfer applications including pharmaceutical, laboratory; bioprocess manufacturing; and food and beverage applications. There are no phthalate plasticizers that can leach out into transported fluid. Manufactured under strict process controls for precise dimensional tolerances. Material complies with Food-Grade, USDA, USP Class VI and 3A standards. Translucent. Non-cytotoxic, non-pyrogenic and non-hemolytic. NOTE: Not intended for invasive use. Autoclavable/USPVI

Cat. No.	I.D. x O.D. x Wall, in.	Max. Press. psig at 73°F; 23°C	Max. Press. psig at 320°F; 160°C	Ft. per Case
8060-0020	1/16 x 1/8 x 1/32	15	9	50
8060-0030	1/8 x 1/4 x 1/16	21	12	50
8060-0040	3/16 x 5/16 x 1/16	13	8	50
8060-0050	3/16 x 3/8 x 3/32	14	8	50
8060-0060	1/4 x 3/8 x 1/16	12	7	50
8060-0070	1/4 x 7/16 x 3/32	17	10	50
8060-0080	1/4 x 1/2 x 1/8	17	10	50
8060-0100	5/16 x 1/2 x 3/32	16	10	50
8060-0120	3/8 x 1/2 x 1/16	11	7	50
8060-0130	3/8 x 9/16 x 3/32	12	7	50
8060-0140	3/8 x 5/8 x 1/8	14	8	50
8060-0190	1/2 x 11/16 x 3/32	10	6	50
8060-0200	1/2 x 3/4 x 1/8	12	7	50
8060-0260	5/8 x 7/8 x 1/8	12	7	50
8060-0310	3/4 x 1 x 1/8	9	5	50

Designed to provide extended pump life in peristaltic applications.

Cat. No.	Fits Pump Tube Size	Hose Barb, in.	Hose Barb, mm	Ft. per Case
8060-3015	#15	3/16	4.8	25
8060-3016	#16	1/8	3.2	25
8060-3017	#17	1/4	6.4	25
8060-3018	#18	5/16	7.9	25
8060-3024	#24	1/4	6.4	25
8060-3025	#25	3/16	4.8	25
8060-3173	#73	3/8	9.5	25
8060-3182	#82	1/2	12.7	25



A

USP VI

8061 NALGENE 65 Braided Platinum-Cured Silicone Tubing

Manufactured with interwoven polyester braid for additional pressure and kink resistance, Shore A 65. Durable, high-purity tubing is designed for a variety of pump and transfer applications including pharmaceutical, laboratory; bioprocess manufacturing; and food and beverage applications. There are no phthalate plasticizers that can leach out into transported fluid. Manufactured under strict process controls for precise dimensional tolerances. Material complies with Food-Grade, USDA, USP Class VI and 3A standards. Translucent. Non-cytotoxic, non-pyrogenic and non-hemolytic. NOTE: Not intended for invasive use. Autoclavable/USPVI

Cat. No.	I.D. x O.D. x Wall, in.	Max. Press. psig at 73°F; 23°C	Max. Press. psig at 320°F; 160°C	Ft. per Case
8061-3070	1/4 x 9/16 x 5/32	156	78	25
8061-3130	3/8 x 21/32 x 5/32	136	68	25
8061-3190	1/2 x 13/16 x 5/32	126	63	25