

PDS No. 2251xx	PRODUCT DATA SHEET				Page 1 of 1	
Revision 03	PS Multipurpose Container, 175 ml				 greiner bio-one	
	Greiner Item-No. 2251xx					
Valid for Item-No.:	225170	225180				

1.	Description / Specification	
1.1	Description	PS Multipurpose Container, 175 ml, with metal screw cap 225170: standard, no label 225180: standard, plain label
1.2	Dimensions	See Customer drawing
1.3	Volume	Working volume: 155 ml
1.4	Material / Resin	Container: PS (Polystyrene), free from heavy metal Screw cap: Tin plate, free from heavy metal Label: Paper, free from heavy metal
1.5	Colour	Container: clear Screw cap: silver 225180: Label: white
1.6	Sterilization	No
1.7	Quality Control	- Raw Material-Control: physical testing - Product-Control: testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	- Single use only - Aseptically produced (Class II Clean Room)

2.	Features	
2.1	Basic features	-
2.2	Temperature range	-20°C to +60°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	Container: see homepage: www.gbo.com/bioscience →Products→Literature →Technical Information →Chemical Resistance of Resins
2.6	Shelf life	N/A
2.7	Other Information	-

3.	Packaging	
3.1	Pieces / Bag	20
3.2	Pieces / Box	120
3.3	Lot-No.	YYWWXX (year, week, day, shift)
3.4	Other Information	-

4.	Other Information	
	-	

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 02	Date 31 July 2014	Date 1 August 2014	Date 1 August 2014	
Date 02.12.2011	Name S. Kaelberer	Name Dr. T. Schreiber	Name A. Schulz	