

Certificate of Quality

VWR Catalogue Number	See List
Description	VWR® Water Sampling Kits
Country of Origin	Manufactured in Canada
Date of Issue (yyyy-mm-dd)	2025-07-30

VWR Catalogue Number	Capacity	Thickness	W x H	Unit
89003-816	180 ml (6 oz.)	2.5 mil.	76 × 178 mm (3 × 7")	Case of 100
89003-818	540 ml (18 oz.)	3 mil.	114 × 229 mm (4.5 × 9")	Case of 100

Quality System Compliance

This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared their Quality Management System (QMS) complies with **ISO 9001:2015** standard.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable product specifications.

Random tests are carried out for aerobic and anaerobic Bacteria as well as yeast and molds. Sterility Tests, based on USP 40 Microbiological Tests / <71> Sterility Tests, are performed to guaranty sterility. of the processes.

Product Specifications

Material	Polyethylene, Polypropylene tape, Steel Wires, Thiosulfate tablets.
Sterilization	Products are sterilized by gamma irradiation and dose released upon ISO11137-1 and 11137-2 recommended practices in effect at the time of validation and meets a minimum requirement of 10 ⁻³ SAL with a specified dose range of 15-32 kGy.
ATP Assay	Not Tested.
DNase & RNase Free	The products are free of any detectable DNase/RNase contamination.

BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not Tested.
Non-Pyrogenic Statement	The products are free from Pyrogen.
Latex Statement	The products are Latex free.
BPA Statement	Bisphenols are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	DEHP is not used in the formulation or manufacture of the other raw material.
RoHS	Not applicable.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list are used in manufacturing of the raw materials used to produce this product. No routine analysis is performed.
Storage Conditions	Store at ambient temperature.
Shelf Life	3 years.

Disclaimer: VWR'S Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability and distribution processes. VWR is not a manufacturer and does not engage in manufacturing, co- manufacturing, or processing activities such as design, packaging, labelling, testing, or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.

VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use. This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoQ using new layout	July 2025