

Certificate of Conformity

VWR Catalogue Number	See List
Description	Sleeves for standard ground joints
Country of Origin	Manufactured in United Kingdom
Date of Issue (yyyy-mm-dd)	2021-05-11

Catalogue Number	For NS
201-0025	10/19
201-0026	12/21
201-0027	14/23
201-0028	19/26
201-0029	24/29
201-0030	29/32
201-0031	34/35
201-0032	40/38
201-0033	45/40
201-0034	50/42
201-0035	55/44
201-0024	7/16

Quality System Compliance

Products are manufactured under the **ISO 9001:2015** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material	Virgin Polytetrafluoroethylene (PTFE) Grade TF1620. Material is in compliance with Directive FDA 21 CFR 177.1550 and complies with USP Class VI - USP 27, NF 22,2004 for plastic at 70 °C.
Sterilization	Non-Sterile.
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Bisphenol are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	Please see RoHS information below.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analysis is performed.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list published 01-17-2021 are used in manufacturing of the raw materials used to produce this product. No routinely analysis is performed.
Storage Conditions	No special storage conditions required.
Shelf Life	Indefinite when stored under ideal conditions.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
 This document has been produced electronically and is valid without a signature.