

Certificate of Conformity

VWR Catalogue Number	See List
Description	Colour-frosted microscope slides
Country of Origin	Manufactured in China
Date of Issue (yyyy-mm-dd)	2025-08-05

Cut edges

Catalogue Number	Designation	Colour
631-1558	Single-frosted slides, coloured ends	Green
631-1559	Single-frosted slides, coloured ends	Orange
631-1555	Single-frosted slides, coloured ends	Blue
631-1557	Single-frosted slides, coloured ends	Yellow
631-1554	Single-frosted slides, coloured ends	White
631-1556	Single-frosted slides, coloured ends	Rose

Bevelled edges

Catalogue Number	Designation	Colour
631-1560	Single-frosted slides, coloured ends	White
631-1562	Single-frosted slides, coloured ends	Blue
631-1565	Single-frosted slides, coloured ends	Orange
631-1561	Single-frosted slides, coloured ends	Rose
631-1563	Single-frosted slides, coloured ends	Yellow
631-1564	Single-frosted slides, coloured ends	Green

Quality System Compliance

This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared their Quality Management System (QMS) complies with **ISO 9001:2015**.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material	Float Glass.
Length	75.50 ~ 76.00 mm
Width	25.50 ~ 26.00 mm
Thickness	1.00 ~ 1.20 mm
Refractive index	1.53 +/-0.02 of refractive index of the light at the green mercury e-line (wavelength=547.07nm) ISO 8037-2.
Wettability	The central parts of the two surfaces of the slides should be wetted equally, and slightly dehydration of the edges is allowed. (ISO 8037-2)
Sterilization	Not applicable.
ATP Assay	Not applicable.
DNase & RNase Free	Not applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	Not applicable.
Latex Statement	Not applicable.
BPA Statement	Not applicable.
DEHP Statement	Not applicable.
RoHS	Not applicable.

REACH Statement	Not applicable.
Storage Conditions	Room temperature, dry.
Shelf Life	No defined shelf life.

Disclaimer: VWR's Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability, and distribution processes. VWR is not a manufacturer and does not engage in manufacturing, co-manufacturing or processing activities such as design, packaging, labeling, testing or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.

VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use.

This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoQ using new layout	04-2022
2	2.0	Modified Quality system compliance and disclaimer	08-2025