

Certificate of Conformity

VWR Catalogue Number	113-8226
Description	Coveralls with Hood, with Elasticated cuffs, Ankles and waist, Thumb Loop, White Colour, Size M, Non-Sterile
Lot Number	201010
Date of Manufacture (yyyy-mm-dd)	2020-10-10
Country of Origin	Manufactured in China
Date of Issue (yyyy-mm-dd)	2022-08-23

Quality System Compliance

Products are manufactured under the **ISO 9001:2015 & ISO 13485:2016** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements. Inspected according to AQL level II 4.0 sampling plan per MIL-STD-105E.

Inspected according to ASTM D5034 standard.

	MD	CD
Break Tensile Strength [N]	152.2	69.1
Elongation rate [%]	124.3	116.5
Trapezoid Tearing Strength (N)	157.2	81.2

Product Specifications

Material	60g Spunbonded Polypropylene fabric/Meltblown /Spunbonded Polypropylene fabric (SMS).
Sterilization	Non-Sterile.
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.

BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Not Applicable.
DEHP Statement	Not Applicable.
RoHS	Not Applicable.
REACH Statement	Not Applicable.
Storage Conditions	Store at 0° to 50°C, 45% to 65%.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.