

Certificate of Conformity

VWR Catalogue Number	See List
Description	VWR®, Septa ND20
Country of Origin	Manufactured in USA - 548-0797A, 548-0807A, 548-1059A, 548-1484A, 548-1979A, 548-1833A, 548-1456A, 548-0570A, 548-0569A, 548-0869A, 548-0572A, 548-0565A, 548-0567A, 548-0568A, 548-0573A, 548-0574A Manufactured in China - 548-0571A
Date of Issue (yyyy-mm-dd)	2023-12-26

Catalogue Number	Septum	Hardness	Thickness [mm]
548-0797A	Red butyl/grey PTFE	50° shore A	1.3
548-0807A	White silicone/blue PTFE	55° shore A	1.5
548-1059A	Transparent blue silicone/transparent PTFE	45° shore A	3.0
548-1484A	Transparent blue silicone/white PTFE	45° shore A	1.3
548-1979A	PTFE virginal	53° shore D	0.2
548-1833A	White silicone/red PTFE	45° shore A	1.3
548-1456A	Transparent blue silicone/white PTFE, Y-slit	45° shore A	1.3
548-0570A	Moulded dark grey butyl	50° shore A	3.0
548-0569A	Moulded grey butyl/PTFE	50° shore A	3.0
548-0869A	Viton 1A black	70° shore A	1.5
548-0572A	White silicone/silver aluminium foil	45° shore A	3.0
548-0565A	Red-orange natural rubber/transparent TEF	60° shore A	1.3
548-0567A	Transparent blue silicone/white PTFE	45° shore A	1.3
548-0568A	Pharma-fix bromo-butyl/PTFE	50° shore A	3.0
548-0571A	Transparent blue silicone/white PTFE	45° shore A	3.0
548-0573A	White silicone/beige PTFE, high temperature quality	45° shore A	3.2
548-0574A	Red-orange natural rubber/transparent TEF	60° shore A	1.3

Quality System Compliance

Products are manufactured under the **ISO 9001:2015** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable product specifications.

Product Specifications

Material	Refer table.
Sterilization	Not applicable.
ATP Assay	Not applicable.
DNase & RNase Free	Not applicable.
BSE/TSE	No use of any raw material produced from, or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	Not applicable.
Latex Statement	Not applicable.
BPA Statement	Not applicable.
DEHP Statement	Not applicable.
RoHS	Not applicable.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list published 2023-06-14 are used in manufacturing of the raw materials used to produce this product. No routinely analyse is performed.

Storage Conditions Not applicable.

Shelf Life Not applicable.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoQ using new layout	12-2023