

Certificate of Conformity

VWR Catalogue Number	See List
Description	Cell slides in culture dishes, Individually Wrapped
Country of Origin	Manufactured in China
Date of Issue (yyyy-mm-dd)	2021-12-21

Catalogue Number	Description	Diameter [mm]	Growth area [cm ²]
734-2973	Culture dish, 60 mm, with 18 cell slides, cell growth area 0,50 cm ² per cell slide	8	9.00
734-2974	Culture dish, 60 mm, with 12 cell slides, cell growth area 0,79 cm ² per cell slide	10	9.42
734-2975	Culture dish, 100 mm, with 45 cell slides, cell growth area 0,50 cm ² per cell slide	8	22.50
734-2976	Culture dish, 100 mm, with 32 cell slides, cell growth area 0,79 cm ² per cell slide	10	25.12

Quality System Compliance

Products are manufactured under the **ISO 9001:2015 & ISO 13485:2016** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

This Product is Class I category device as defined by the FDA in 21 CFR Parts 862-892.

Material Polystyrene (PS) + Polyethylene terephthalate (PET).

Food contact statement The composition of the product complies with the requirements of the 21 CFR 177.1640 "Polystyrene and rubber modified polystyrene".

Sterilization	Product labelled as sterile is EB (Electron beam) irradiated as per ISO 11137 and meets SAL 10 ⁻⁶ with a specified dose range of 15-30 kGy.
ATP Assay	Not Applicable.
DNase & RNase Free	This product is free of any detectable DNase/RNase contamination.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Testing is conducted to quality all material resins using ISO 10993 standards for cytotoxicity and have been shown to be non-toxic.
Non-Pyrogenic Statement	The acceptance level for product is 0.5 EU/mL or less than 20 EU/device (TAL Gel Clot Method).
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Bisphenol are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	Please see RoHS information below.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.
REACH Statement	Not Applicable.
Storage Conditions	Store at normal temperature.
Shelf Life	3 years.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
 This document has been produced electronically and is valid without a signature.