



## ZPACK13951000 PLASTIC DRUM LIQ 10 L DIN 51 UN



### SPECIFICATIONS

|                    |               |
|--------------------|---------------|
| UN-certification : | 3H1/Y1.9/200/ |
| Material :         | HDPE          |
| Neck Type :        | DIN 51        |
| Height :           | 310,00 mm     |
| Length :           | 230,00 mm     |
| Wide :             | 198,00 mm     |

### AVAILABLE ACCESSORIES

|                   |           |
|-------------------|-----------|
| Key for opening : | 29548.064 |
| Tap DIN 51 :      | 29548.042 |

Remark : Packaging materials and specifications shown may be modified without preceding notice  
Version dd. 22 mai 2012

tel/tél +32(0)16 385 011  
fax +32(0)16 385 385  
[info@be.vwr.com](mailto:info@be.vwr.com)  
[www.vwr.com](http://www.vwr.com)

Researchpark Haasrode 2020  
Geldenaaksebaan 464 B-3001 Leuven  
RPR/RPM Leuven  
BTW/TVA BE 0403.593.343