

Certificate of Conformity

VWR Catalogue Number	See List
Description	VWR®, Safety Waste Caps, with tube connections
Country of Origin	Manufactured in Germany
Date of Issue (yyyy-mm-dd)	2023-11-09

Catalogue Number	Description	Thread	Connector
590-0043	Safety waste cap, including connector for exhaust filter	GL 45	3× for tubing 1.6/3.2 mm OD
590-0044	Safety waste cap, including connector for exhaust filter	GL 45	2× for tubing 1.6/3.2 mm OD, 1× for tubing 6.4 - 9 mm OD
590-0045	Safety waste cap, including connector for exhaust filter	S 55	2× for tubing 1.6/3.2 mm OD, 1× for tubing 6.4 - 9 mm OD
590-0046	Safety waste cap, including connector for exhaust filter	S 60/61	2× for tubing 1.6/3.2 mm OD, 1× for tubing 6.4 - 9 mm OD

Quality System Compliance

Products are manufactured under the **DIN EN ISO 9001:2015** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material	PTFE (Polytetrafluoroethylene) PE (Polyethylene)
Sterilization	Non-Sterile.
ATP Assay	Not Applicable.

DNase & RNase Free	This product is free of any detectable DNase/RNase contamination.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Bisphenol are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	Not Applicable.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list published 2023-06-14 are used in manufacturing of the raw materials used to produce this product. No routinely analyse is performed.
Storage Conditions	Room temperature.
Shelf Life	Unlimited.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoQ using new layout	11-2023