

Sodium Hydroxide
Baker Analyzed
Pellets



Material No: 0348.5000
Effective from: 2019/06/05

Product Specification

Meets ACS Specifications.

Test	Specification
Assay	$\geq 97.0 \%$
Calcium (Ca)	$\leq 0.005 \%$
Chloride (Cl)	$\leq 0.006 \%$
Heavy Metals (as Ag)	$\leq 0.001 \%$
Magnesium (Mg)	$\leq 0.002 \%$
Potassium (K)	$\leq 0.02 \%$
Sodium Carbonate (Na_2CO_3)	$\leq 1.0 \%$
Trace Impurities (in ppm)	
Iron (Fe)	≤ 3
Mercury (Hg)	≤ 0.1
Nickel (Ni)	≤ 5
Nitrogen Compounds (as N)	≤ 3
Phosphate (PO_4)	≤ 2
Sulphate (SO_4)	$\leq 30 \text{ ppm}$

For Laboratory, Research or Manufacturing use.

Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000