

Material	28245.265
Material description	Sodium hydroxide, pellets
Grade	GPR RECTAPUR
CAS Number	1310-73-2
Molecular formula	NaOH
Molecular mass	40.00

Characteristics	Specifications
Assay	≥ 98 %
Heavy metals (as Pb)	≤ 5 ppm
Total N (Nitrogen)	≤ 5 ppm
CO ₃ (as Na ₂ CO ₃)	≤ 1.0 %
Cl (Chloride)	≤ 0.02 %
PO ₄ (Phosphate)	≤ 20 ppm
SiO ₄ (as SiO ₂)	≤ 100 ppm
SO ₄ (Sulphate)	≤ 50 ppm
Al (Aluminium)	≤ 10 ppm
Fe (Iron)	≤ 10 ppm
K (Potassium)	≤ 0.1 %
Pb (Lead)	≤ 5 ppm

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Anja Vanhalle, Head of Laboratory - Haasrode
VWR International bv; Geldenaaksebaan 464; BE-3001
Leuven; Belgium

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.