

Certificate of Conformity

VWR Catalogue Number	See List
Description	Mortars without Pestles, with Spout
Country of Origin	Manufactured in Czech Republic
Date of Issue (yyyy-mm-dd)	2021-10-25

Catalogue Number	Capacity [ml]	Height [mm]	Internal Diameter [mm]
410-0108	20	32	50
410-0109	70	42	63
410-0110	75	46	76
410-0111	140	50	85
410-0112	175	56	92
410-0113	250	64	115
410-0114	400	71	122
410-0115	650	90	155
410-0116	1000	104	185
410-0117	2400	127	255
410-0118	4300	160	265

Quality System Compliance

Products are manufactured under the **ISO 9001:2015** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements. Inspected according to AQL level 2.5 sampling plan. Thermal shock resistance as per ASTM & Dimension accuracy.

Product Specifications

Material	Porcelain mass C110. The material grades conform to EU Directive 1935/2004 on materials and articles intended to come into contact with food.
Sterilization	Not applicable.
ATP Assay	Not applicable.
DNase & RNase Free	Not applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/ (TSE).
Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	Not applicable.
Latex Statement	Not applicable.
BPA Statement	Not applicable.
DEHP Statement	Not applicable.
RoHS	Not applicable.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list published 2021-07-08 are used in manufacturing of the raw materials used to produce this product. No routinely analysis is performed.
Storage Conditions	No Special Storage Condition Required.
Shelf Life	Unlimited.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.