




Technical Datasheet

VWR® FACE MASK 3-PLY EARLOOPS

CHARACTERISTICS:

- Made of non-woven fabric 20g PP/20g MB/25g PP
- Multiple colours option from white, blue and green.
- Ear loops
- Bacteria filtration rate >90%
- Nature Rubber Latex Free
- Comply with European Commission Directive 1935/2004/EC, Regulation (EU) No.10/2011, Regulation (EU) No.2016/1416 and Regulation (EU) 2017/752.

PRODUCTS AND DIMENSION SPECIFICATION



Unit:CM

ECN	Colour	Dimension			
		A	B	C	D
113-8927	White	17.5	9.5	12	17.5
113-0716	Blue	17.5	9.5	12	17.5
113-0714	Green	17.5	9.5	12	17.5
Tolerance		±1	±1	±1	±1

Technical Datasheet

PACKING:

Package/Quantity	50pcs/box, 20boxes/case
Best Storage Conditions	Temp : 0- 50°C, Humidity : 45-85% without sunshine exposure
Country of Origin	China

QUALITY ASSURANCE CRITERIA:

1. Dimension and Function - Critical (AQL1.5, S-2)

- Product dimension out of spec.
- Fabric characteristics out of spec

2. Visual Inspection/ Attributes - Major (AQL 2.5, L-II)

- Foreign Material
- Dirt
- Thread ends loose or too long(over 1cm)
- Holes/Tears
- Foul odor
- Discolor
- Cloth surface applied wrongly

3. Visual Inspection/Attributes - Minor (AQL 4.0, L-II)

- Unwanted folds
- Elastic tightness non-consistent

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **November 2015**