



Specification

1.00045.0010 MUP Selective Supplement

Specification

Appearance (clearness) clear

Testing as supplement in TOS-Propionat Agar (1.00043.0500).

Specification

Recovery on test medium (Bifidobacterium animalis subsp. lactis ATCC 27536)	≥ 70	%
Recovery on test medium (Bifidobacterium breve ATCC 15700)	≥ 70	%
Recovery on test medium (Bifidobacterium longum ATCC 15707)	≥ 70	%
Growth (Lactobacillus casei ATCC 393 (WDCM 00100))	no growth	
Growth (Lactobacillus delbrueckii subsp. bulgaricus ATCC 11842 (WDCM 00102))	no growth	
Growth (Streptococcus salivarius subsp. thermophilus DSM 20259)	no growth	

Incubation: 72 hrs. at 37°C, anaerobic.

Stefanie Fischer

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.