

Material	24975.294
Material description	D(-)-Fructose
Grade	GPR RECTAPUR
CAS Number	57-48-7
Molecular formula	C ₆ H ₁₂ O ₆
Molecular mass	180.16

Characteristics	Specifications
Specific optical rotation (10 %; water)	-94 - -90 °
Ignition residue (SO ₄)	≤ 0.1 %
Water	≤ 0.5 %
Cl (Chloride)	≤ 100 ppm
SO ₄ (Sulphate)	≤ 100 ppm

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Jessie Pijck, Head of Laboratory - Haasrode
VWR International bv; Geldenaaksebaan 464; BE-3001
Leuven; Belgium

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.