

Product Specification Sheet

VWR CELL SCRAPERS

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

Cat. No.	Description	Blade length	Total length	Pk
734-2602	Cell scraper,scraper blade	20 mm	250 mm	100
734-2603	Cell scraper,lifter blade	20 mm	250 mm	100
734-2604	Cell scraper,scraper blade	30 mm	390 mm	100
734-2605	Cell scraper,lifter blade	30 mm	390 mm	100

Features:

- Temperature range: -20°C to +60°C
- Autoclavability: No
- Non-pyrogenic
- Shelf life 3 years


Sterility: Sterilized by gamma irradiation SAL10⁻⁶

Other Characteristics:


- Specifically designed to make collecting of cells easier and more effective
- Cross-ribbed handle in ABS provides a greater rigidity and ensures a better control while scraping cells
- Thin, flexible 20 mm blade in TPE prevents cells from any damage
- Individually wrapped
- Manufactured in a class 100,000 room environment, under EN ISO13485:2003 and ISO9001:2008 quality management system

Product Specification Sheet

VWRI734-2602 and VWRI734-2603:



VWRI734-2604 and VWRI734-2605:



Disclaimer: VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there is any doubt, ask us.

Date of issue: **March 2014**