



ZPACK13962500 PLASTIC BOTTLE LIQ 2,5 L DIN 45 UN X



SPECIFICATIONS

UN-certification : 3H1/X1.3/Y1.9/250/..... Material : HDPE

Neck Type : DIN 45

Height : 293,00 mm

Length : 134,00 mm

Wide : 134,00 mm

AVAILABLE ACCESSORIES

Adaptator for dispenser S40 to GL 45 : 217-9271

Key for opening : 84544.001

Remark : Packaging materials and specifications shown may be modified without preceding notice
Version dd. 22 mai 2012

tel/tél +32(0)16 385 011
fax +32(0)16 385 385
info@be.vwr.com
www.vwr.com

Researchpark Haasrode 2020
Geldenaaksebaan 464 B-3001 Leuven
RPR/RPM Leuven
BTW/TVA BE 0403.593.343