



TECHNICAL DATA SHEET

Ideal for shipping, storing, or staining slides. These polypropylene containers hold five standard slides (75 x 25 mm and 3 x 1"). Inside wall features slotted channels to keep slides separated and to make removal easier. Close securely to keep slides safe. The end-opening mailer is 43 x 24 x 88 mm. Pack of 25.

Product name:	5-Place Slide Mailer, PP, end opening, Natural pk 25
Part number:	HS15986
Color:	Natural
Material of construction	Polypropylene
Unit Dimensions	43 x 24 x 88 mm
Packaging:	6.8x4x4.3" (inches)

Heathrow Scientific®, LLC

620 Lakeview Parkway ▪ Vernon Hills ▪ Illinois ▪ 60061 ▪ Phone: (847) 816-5070 ▪ Fax: (847) 816-5072
E-mail: info@heathrowscientific.com ▪ Website: www.heathrowscientific.com