



Certificate of Conformity

VWR DIGITAL TIMER, ONE CHANNEL, TRACEABLE® CLIP-IT

VWR International certifies that the products comply in every respect with the relevant particular specifications and/or drawings and/or relevant standards and regulations in force, mentioned in this document.

VWR International certifie que les produits sont conformes à tous égards aux spécifications particulières appropriées et / ou dessins, normes et règlements en vigueur, mentionnées dans le présent document

VWR International bestätigt, dass das Produkte allen zutreffenden Spezifikations- und/oder Zeichnungs- und/oder Standardanforderungen sowie den im Dokument aufgeführten gültigen Vorschriften entspricht.

Cat No: **609-0205**

Accuracy: **±0.005 %**

Size: **50 x 19 x 38 mm**

Resolution: **1 second**

Colour: **White/Blue**

Other Characteristics:

- Count up and down.
- Capacity of 100 hours in steps of one second.
- At zero, an alarm sounds for 60 seconds or may be silenced manually.
- Helpful memory saves last programmed setting for recall in repetitive events.
- Complies with ISO/IEC 17025, ANSI/NCSL Z540-1 and 9001, accredited by A2LA indicates traceability to standards provided by NIST.

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **February 2017**