

Certificate of Quality

VWR Catalogue Number	See List
Description	VWR® Self-Seal Sterilization Pouches
Country of Origin	Manufactured in Canada / USA
Date of Issue (yyyy-mm-dd)	2025-07-22

VWR Catalogue Number	Dimensions
89140-800	13.3 × 25.4 cm (5¼ × 10")
89140-802	19.1 × 33 cm (7½ × 13")
89140-804	30.5 × 45.7 cm (12 × 18")

Quality System Compliance

This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared their Quality Management System (QMS) complies with **ISO 13485:2016** standard.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable product specifications.

Product Specifications

Material	Pouch paper Medical grade paper 60 GSM Laminate film CPP/48ga, PET Adhesive tape Made of polyester film coated with medical-grade acrylic adhesive.
Sterilization	Not applicable.
ATP Assay	Not applicable.
DNase & RNase Free	The products are free of any detectable DNase/RNase contamination.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).

Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	Not applicable.
Latex Statement	The products are not made with natural rubber latex.
BPA Statement	Bisphenols are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	Not applicable.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP) are used in manufacturing the raw materials and final product. No routine analysis is performed.
REACH Statement	Not applicable.
Storage Conditions	Store at room temperature of 15 °C (59 °F) and 25 °C (77 °F) with relative humidity of 40% - 70%. Storage location should be clean, dust free, and away from direct sunlight, fluorescent or ultraviolet light.
Shelf Life	3 years.

Disclaimer: VWR'S Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability and distribution processes. VWR is not a manufacturer and does not engage in manufacturing, co- manufacturing, or processing activities such as design, packaging, labelling, testing, or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.

VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use. This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoQ using new layout	September 2021
2	1.0	Updated Cytotoxicity statement (Internal error)	October 2024
3	1.0	Updated incorrect dimensions	February 2025
4	1.1	Updated RoHS statement and Country of origin	May 2025
5	2.0	Updated ISO and Disclaimer Statements	July 2025