

Liquid handling & sample preparation consumables



Table of contents

01

Pipetting Solution
Pages 3-4

02

Tips
Pages 5-8

03

Tubes
Pages 9-10

04

PCR Consumables
Pages 10-11

05

Reservoirs
Pages 11

06

Loops & T-Spreader
Pages 11

07

Dispenser
Pages 12

PHP (PRECISE HIGH PERFORMANCE) PIPETTE

SINGLE & MULTI CHANNEL PIPETTES, MECHANICAL, PRECISE HIGH PERFORMANCE (PHP) MICRO PIPETTE

Key features:

- Magnet assisted piston for consistently precise results
- Innovative spring and seal design for one of the lowest plunge force requirement
- Durable tip ejector made using special corrosion resistant alloy, protecting against chemical & physical corrosion
- Effortless tip ejection with unique shock absorbing mechanism
- Spring loaded tip cone for lesser strain & proper tip fitment for precise pipetting (for Multi-Channel only)
- Fully Autoclavable
- Ergonomic design with proper weight distribution
- Calibrated as per ISO 8655-6 in ISO 17025 accredited lab
- Large and clear 4 digit display
- Color coding for easy identification
- Easy cleaning & in-house calibration
- Highly durable PVDF tip cone
- 3-Year warranty

Ordering information:

VWR PHP Pipettes

Description	UOM	Cat No.	Price
PHP - SINGLE CHANNEL VV VV PIPETTE 0.1-2.5	1*1 item	VWRA613-6688	14500
PHP - SINGLE CHANNEL VV PIPETTE 0.5-10	1*1 item	VWRA613-6689	14500
PHP - SINGLE CHANNEL VV PIPETTE 2-20	1*1 item	VWRA613-6690	14500
PHP - SINGLE CHANNEL VV PIPETTE 5-50	1*1 item	VWRA613-6691	14500
PHP - SINGLE CHANNEL VV PIPETTE 10-100	1*1 item	VWRA613-6692	14500
PHP - SINGLE CHANNEL VV PIPETTE 20-200	1*1 item	VWRA613-6693	14500
PHP - SINGLE CHANNEL VV PIPETTE 100-1000	1*1 item	VWRA613-6694	14500
PHP - SINGLE CHANNEL VV PIPETTE 1000-5000	1*1 item	VWRA613-6695	14500
PHP - SINGLE CHANNEL VV PIPETTE 1000-10000	1*1 item	VWRA613-6696	14500
PHP - SINGLE CHANNEL FIX VOL PIPETTE 2.5	1*1 item	VWRA613-6697	10000
PHP - SINGLE CHANNEL FIX VOL PIPETTE 5	1*1 item	VWRA613-6698	10000
PHP - SINGLE CHANNEL FIX VOL PIPETTE 10	1*1 item	VWRA613-6699	10000
PHP - SINGLE CHANNEL FIX VOL PIPETTE 20	1*1 item	VWRA613-6700	10000
PHP - SINGLE CHANNEL FIX VOL PIPETTE 25	1*1 item	VWRA613-6701	10000
PHP - SINGLE CHANNEL FIX VOL PIPETTE 50	1*1 item	VWRA613-6702	10000
PHP - SINGLE CHANNEL FIX VOL PIPETTE 100	1*1 item	VWRA613-6703	10000
PHP - SINGLE CHANNEL FIX VOL PIPETTE 200	1*1 item	VWRA613-6704	10000
PHP - SINGLE CHANNEL FIX VOL PIPETTE 250	1*1 item	VWRA613-6705	10000
PHP - SINGLE CHANNEL FIX VOL PIPETTE 500	1*1 item	VWRA613-6706	10000

Description	UOM	Cat No.	Price
PHP - SINGLE CHANNEL FIX VOL PIPETTE 1000	1*1 item	VWRA613-6707	10000
PHP - SINGLE CHANNEL FIX VOL PIPETTE 2000	1*1 item	VWRA613-6708	12500
PHP - SINGLE CHANNEL FIX VOL PIPETTE 5000	1*1 item	VWRA613-6709	12500
PHP - SINGLE CHANNEL FIX VOL PIPETTE 10000	1*1 item	VWRA613-6710	12500
PHP 8 Channel Pipette 0.5 - 10 µl	1*1 item	VWRA613-6711	27000
PHP 8 Channel Pipette 2 - 20 µl	1*1 item	VWRA613-6712	27000
PHP 8 Channel Pipette 5- 50 µl	1*1 item	VWRA613-6713	27000
PHP 8 Channel Pipette 10 - 100 µl	1*1 item	VWRA613-6714	27000
PHP 8 Channel Pipette 20 - 200 µl	1*1 item	VWRA613-6715	27000
PHP 8 Channel Pipette 30 - 300 µl	1*1 item	VWRA613-6716	27000
PHP 12 Channel Pipette 0.5 - 10 µl	1*1 item	VWRA613-6717	32000
PHP 12 Channel Pipette 2 - 20 µl	1*1 item	VWRA613-6718	32000
PHP 12 Channel Pipette 5 - 50 µl	1*1 item	VWRA613-6719	32000
PHP 12 Channel Pipette 10 - 100 µl	1*1 item	VWRA613-6720	32000
PHP 12 Channel Pipette 20 - 200 µl	1*1 item	VWRA613-6721	32000
PHP 12 Channel Pipette 30 - 300 µl	1*1 item	VWRA613-6722	32000
PHP Carousel stand for 6 pipettes	1*1 item	VWRA613-6723	3800
PHP Shelf mounting stand for Pipette	1*1 item	VWRA613-6724	550



A- SERIES PIPETTES

Fully autoclavable Click stop counter system for convenient volume setting; Ergonomic design ensuring light weight & soft plunger movement.

- Durable tip cone made of PVDF offers the highest resistance to physical shocks as well as to chemical corrosion and is compatible with most types of tips
- Calibrated as per ISO 8655-6 in ISO 17025 accredited lab
- 1-Year warranty

Ordering information:

Description	UOM	Cat No.	Price
Fully Auto Var Vol Micropipette 0.5-10µl	1*1 items	VWRA613-6591	3927
Fully Auto Var Vol Micropipette 5-50µl	1*1 items	VWRA613-6592	3927
Fully Auto Var Vol Micropipette 50-200µl	1*1 items	VWRA613-6593	3927
Fully Auto Var Vol pipette 200-1000µl	1*1 items	VWRA613-6594	3927
Fully Auto Var Vol Micropipette 10-100µl	1*1 items	VWRA613-6595	3927
Fully Auto Var Vol pipette 100-1000µl	1*1 items	VWRA613-6596	3927
Fully Auto Var Vol Micropipette 2-20µl	1*1 items	VWRA613-6597	3927
Fully Auto Var Vol Micropipette 20-200µl	1*1 items	VWRA613-6598	3927
Fully Auto Var Vol pipette 1000-5000µl	1*1 items	VWRA613-6599	4488
Fully Auto Var Vol pipette 500-5000µl	1*1 items	VWRA613-6600	4488
Fully Auto Var Vol pipette 1000-10000µl	1*1 items	VWRA613-6601	4488
Fully Auto Var Vol pipette 500-10000µl	1*1 items	VWRA613-6602	4488
Fully Auto Fix Volume Micropipette 5µl	1*1 items	VWRA613-6603	2020
Fully Auto Fix Volume Micropipette 10µl	1*1 items	VWRA613-6604	2020
Fully Auto Fix Volume Micropipette 20µl	1*1 items	VWRA613-6605	2020
Fully Auto Fix Volume Micropipette 25µl	1*1 items	VWRA613-6606	2020
Fully Auto Fix Volume Micropipette 50µl	1*1 items	VWRA613-6607	2020
Fully Auto Fix Volume Micropipette 100µl	1*1 items	VWRA613-6608	2020



Description	UOM	Cat No.	Price
Fully Auto Fix Volume Micropipette 200µl	1*1 items	VWRA613-6609	2020
Fully Auto Fix Volume Micropipette 250µl	1*1 items	VWRA613-6610	2581
Fully Auto Fix Volume Micropipette 500µl	1*1 items	VWRA613-6611	2581
Fully Auto Fix Vol Micropipette 1000µl	1*1 items	VWRA613-6612	2581
Fully Auto Fix Vol Micropipette 2000µl	1*1 items	VWRA613-6613	3500
Fully Auto Fix Vol Micropipette 5000µl	1*1 items	VWRA613-6614	3500
Fully Auto Fix Vol Micropipette 10000µl	1*1 items	VWRA613-6615	3500
Fully Auto 8 Ch Micropipette 5 - 50 µl	1*1 items	VWRA613-6616	19000
Fully Auto 8 Ch Micropipette 10 - 100 µl	1*1 items	VWRA613-6617	19000
Fully Auto 8 Ch Micropipette 20 - 200 µl	1*1 items	VWRA613-6618	19000
Fully Auto 8 Ch Micropipette 30 - 300 µl	1*1 items	VWRA613-6619	19000
Fully Auto 8 Ch Micropipette 0.5 - 10 µl	1*1 items	VWRA613-6620	19000
Fully Auto 12 Ch Micropipette 5 - 50 µl	1*1 items	VWRA613-6621	24500
Fully Auto 12 Ch pipette 10 - 100 µl	1*1 items	VWRA613-6622	24500
Fully Auto 12 Ch pipette 20 - 200 µl	1*1 items	VWRA613-6623	24500
Fully Auto 12 Ch pipette 30 - 300 µl	1*1 items	VWRA613-6624	24500
Fully Auto 12 Ch pipette 0.5 - 10 µl	1*1 items	VWRA613-6625	24500

ELECTRONIC PIPETTE – CLASSIC

Key features:

Intuitive graphical interface makes it operationally easy to use

- Lightweight construction makes pipetting effortless and pain-free even after extensive usage
- Auto-sensing magnetic piston for precise pipetting; State-of-the-art stepper motor for long lasting performance
- 5 different pipetting modes for all pipetting needs in a lab: Forward pipetting, Reverse pipetting, Multi dispensing, Mixing, Diluting
- Calibrated as per ISO 8655-6 standards in ISO 17025 accredited lab; 3 point auto-calibration provides most accurate and precise delivery of liquid throughout the volume range
- Ultra-efficient battery allows 5 hours of continuous usage after every full recharge; Intelligent sleep mode conserves battery when not in use
- A smart charging stand which can be used for charging or to rest the pipette when not in use
- Min Vol - 0.5 µl & Max Vol - 5,000 µl



Description	UOM	Cat No.	Price
Electronic Pipette Single Ch 0.5-10µl	1*1 item	VWRA613-6671	32000
Electronic Pipette Single Ch 5-100µl	1*1 item	VWRA613-6672	32000
Electronic Pipette Single Ch 10-200µl	1*1 item	VWRA613-6673	32000
Electronic Pipette Single Ch 50-1000µl	1*1 item	VWRA613-6674	32000
Electronic Pipette Single Ch 250-5000µl	1*1 item	VWRA613-6675	32000
Electronic Pipette 8 Channel 0.5-10	1*1 item	VWRA613-6676	48000
Electronic Pipette 8 Channel 2-20	1*1 item	VWRA613-6677	48000
Electronic Pipette 8 Channel 5-100	1*1 item	VWRA613-6678	48000

Description	UOM	Cat No.	Price
Electronic Pipette 8 Channel 10-200	1*1 item	VWRA613-6679	48000
Electronic Pipette 8 Channel 15-300	1*1 item	VWRA613-6680	48000
Electronic Pipette 12 Channel 0.5-10	1*1 item	VWRA613-6681	52500
Electronic Pipette 12 Channel 2-20	1*1 item	VWRA613-6682	52500
Electronic Pipette 12 Channel 5-100	1*1 item	VWRA613-6683	52500
Electronic Pipette 12 Channel 10-200	1*1 item	VWRA613-6684	52500
Electronic Pipette 12 Channel 15-300	1*1 item	VWRA613-6685	52500
Electronic Pipette Charging Stand	1*1 item	VWRA613-6686	5600
Acrylic Stand for Electronic Pipette	1*1 item	VWRA613-6687	2500

DNase, RNase ENDO TOXIN FREE PRECISE SHAPED TIPS TO SUIT MOST SENSITIVE PIPETTING APPLICATIONS IN LAB

- Highly transparent tips, produced through high quality diamond finish mould and special FDA approved medical grade virgin Polypropylene
- Innovative design for compatibility with a range of pipettes, ensuring minimum force requirement during loading and ejection
- Ideal match for high sensitive PCR works
- Produced in a fully automatic, human touch free, clean room production facility
- Free from RNase, DNase, Human DNA, Pyrogen & PCR inhibitor. Free of Lubricants and heavy metals
- Graduation marks on all the tips that act as visual checks to minimize any risk of errors during sample aspiration
- All these high quality tips are also available as sterile



BULK TIPS

- Packed in easily tearable Polyethylene self standing bag for ease of handling
- Double sealing through perforation and zip-lock for better protection against contamination and spillage



Ordering information:

Bulk Pack:

Description	UOM	Cat No.	Price
10µl uni Grad Bulk tip	1*1000 items	VWRA612-6332	561
10µl Long uni Grad tip	1*1000 items	VWRA612-6388	561
200µl uni Grad tip	1*1000 items	VWRA612-6478	561
200µl Long uni Grad tip, 300ul	1*1000 items	VWRA612-6509	673
1000µl uni Grad tip	1*1000 items	VWRA612-6526	738
1250µl uni Grad tip	1*500 items	VWRA612-6577	393

Sterile Bulk Pack:

Description	UOM	Cat No.	Price
10µl uni Grad tip, Sterile	1*1000 items	VWRA612-6333	833
10µl Long uni Grad tip, Sterile	1*1000 items	VWRA612-6389	833
10µl Narrow uni Grad tip	1*1000 items	VWRA612-6409	561
10µl Narrow uni Grad tip, 20µl . Sterile	1*1000 items	VWRA612-6414	833
200µl uni Grad tip, Sterile	1*1000 items	VWRA612-6479	833
200µl Long uni Gr tip, 300ul ., Sterile	1*1000 items	VWRA612-6511	914
1000µl uni Grad tip, Sterile	1*1000 items	VWRA612-6529	974
1250µl uni Grad tip, Sterile	1*500 items	VWRA612-6213	577

RACKED TIPS

- Fully autoclavable & stackable
- Hinge rack for single hand operation
- Sturdy and solid design for long durability
- Color coded tray for easy identification



Ordering information:

Racked Tips:

Description	UOM	Cat No.	Price
10µl uni Grad tip	1*960 items	VWRA612-6354	2039
10µl Long uni Grad tip	1*960 items	VWRA612-6403	2039
10µl Narrow uni Grad tip	1*960 items	VWRA612-6424	2039
200µl uni Grad tip	1*960 items	VWRA612-6486	2039
200µl Long uni Grad tip, 300ul .	1*960 items	VWRA612-6514	2244
1000µl uni Grad tip	1*576 items	VWRA612-6534	1387
1250µl uni Grad tip	1*576 items	VWRA612-6579	1600

Sterile Racked Tips

Description	UOM	Cat No.	Price
10µl uni Grad tip, Sterile	1*960 items	VWRA612-6370	2250
10µl Long uni Grad tip, Sterile	1*960 items	VWRA612-6404	2250
10µl Narrow uni Grad tip, 20µl . Sterile	1*960 items	VWRA612-6432	2250
200µl uni Grad tip, Sterile	1*960 items	VWRA612-6488	2350
200µl Long uni Gr tip, 300ul ., Sterile	1*960 items	VWRA612-6520	2400
1000µl uni Grad tip, Sterile	1*576 items	VWRA612-6536	1500
1250µl uni Grad tip, Sterile	1*576 items	VWRA612-6220	1803

STERILE LOW RETENTION TIPS

Key fetures:

Specially developed resins combined with high quality diamond finish mould creates super hydrophobic internal surface which ensures maximum recovery of sample yield

- Minimum retention = Maximum accuracy
- No Silicon additives = No leaching

Ordering information:

Sterile Low-retention Bulk

Description	UOM	Cat No.	Price
10µl uni Grad tip, Sterile, LR	1*1000 items	VWRA612-6339	1043
10µl Long uni Grad tip, Sterile, LR	1*1000 items	VWRA612-6122	1043
10µl Narrow uni Gr tip, 20µl . st, LR	1*1000 items	VWRA612-6128	1043
200µl uni Grad tip, Sterile, LR	1*1000 items	VWRA612-6483	1043
200µl Long uni Gr tip, 300ul ., st, LR	1*1000 items	VWRA612-6184	1233
1000µl uni Grad tip, Sterile, LR	1*1000 items	VWRA612-6194	1475
1250µl uni Grad tip, Sterile, LR	1*500 items	VWRA612-6282	1459

Sterile Low-retention Racked

Description	UOM	Cat No.	Price
10µl uni Grad tip, Sterile, LR	1*960 items	VWRA612-6372	2569
10µl Long uni Grad tip, Sterile, LR	1*960 items	VWRA612-6124	3000
10µl Narrow uni Gr tip, 20µl . st, LR	1*960 items	VWRA612-6137	3000
200µl uni Grad tip, Sterile, LR	1*960 items	VWRA612-6491	3100
200µl Long uni Gr tip, 300ul ., st, LR	1*960 items	VWRA612-6190	3100
1000µl uni Grad tip, Sterile, LR	1*576 items	VWRA612-6196	2000
1250µl uni Grad tip, Sterile, LR	1*576 items	VWRA612-6291	2200



FILTER BULK TIPS

- Specially developed PE filters for pipetting application
- The filter serves as a barrier against aerosols formed during pipetting
- Filter tips protect the pipette shafts from contamination and reduce the risk of cross-contamination and faulty measuring results


Ordering information:

Filter Bulk

Description	UOM	Cat No.	Price
10µl uni Grad tip, Filter	1*1000 items	VWRA612-6334	2450
10µl Long uni Grad tip, Filter	1*1000 items	VWRA612-6399	2459
10µl Narrow uni Grad tip, 20µl . Filter	1*1000 items	VWRA612-6419	2459
20µl uni Grad tip, Filter	1*1000 items	VWRA612-6458	2459
50µl uni Grad tip, Filter	1*1000 items	VWRA612-6464	2459
100µl uni Grad tip, Filter	1*1000 items	VWRA612-6470	2459
200µl uni Grad tip, Filter	1*1000 items	VWRA612-6481	2459
200µl Long uni Grad tip, 300ul .,Filter	1*1000 items	VWRA612-6513	2492
1000µl uni Grad tip, Filter	1*1000 items	VWRA612-6531	2951
1250µl uni Grad tip, Filter	1*500 items	VWRA612-6214	1750

Sterile Filter Bulk

Description	UOM	Cat No.	Price
10µl uni Grad tip, Sterile, Filter	1*1000 items	VWRA612-6348	2625
10µl Long uni Grad tip Sterile, Filter	1*1000 items	VWRA612-6123	2625
10µl Narrow uni Gr tip, 20µl . st, flt	1*1000 items	VWRA612-6129	2625
20µl uni Grad tip, Sterile, Filter	1*1000 items	VWRA612-6459	2650
50µl uni Grad tip, Sterile, Filter	1*1000 items	VWRA612-6465	2650
100µl uni Grad tip, Sterile, Filter	1*1000 items	VWRA612-6471	2650
200µl uni Grad tip, Sterile, Filter	1*1000 items	VWRA612-6484	2650
200µl Long uni Gr tip, 300ul ., st, flt	1*1000 items	VWRA612-6189	2650
1000µl uni Grad tip, Sterile, Filter	1*1000 items	VWRA612-6195	3200
1250µl uni Grad tip, Sterile, Filter	1*500 items	VWRA612-6283	1900



Buying on VWR Lab Supplies is as easy as 1-2-3-4



1. Register



[Register](#)
[Forgot Password](#)

2. Explore

Browse Categories

A B C D E F G H I J

Adhesives

Aluminum Foil

Ampules

3. Select

VWR® PHP pipettes

VWR Catalog Number VWRA613-6688

Unit 1*1 items

Price ₹14,500.00

Quantity

4. Order

My Cart

Standard Products

VWR® PHP pipettes

VWRA613-6688

1*1 items

₹14,500.00

VWR, part of Avantor

₹14,500.00

FILTER RACKED TIPS

Key features:

Made of US FDA approved virgin Poly Propylene. Especially designed PE filter to prevent cross contamination. Rack tips support hands-free tip fitting.

Ordering information:

Filter racked Tips

Description	UOM	Cat No.	Price
10µl uni Grad tip, Filter	1*960 items	VWRA612-6371	3300
10µl Long uni Grad tip, Filter	1*960 items	VWRA612-6405	3300
10µl Narrow uni Grad tip, 20µl . Filter	1*960 items	VWRA612-6433	3300
20µl uni Grad tip, Filter	1*960 items	VWRA612-6460	3300
50µl uni Grad tip, Filter	1*960 items	VWRA612-6466	3300
100µl uni Grad tip, Filter	1*960 items	VWRA612-6472	3300
200µl uni Grad tip, Filter	1*960 items	VWRA612-6489	3300
200µl Long uni Grad tip, 300ul ., Filter	1*960 items	VWRA612-6521	3675
1000µl uni Grad tip, Filter	1*576 items	VWRA612-6539	2462
1250µl uni Grad tip, Filter	1*576 items	VWRA612-6221	2624



Sterile Filter Racked

Description	UOM	Cat No.	Price
10µl uni Grad tip, Sterile, Filter	1*960 items	VWRA612-6373	3600
10µl Long uni Grad tip, Sterile, Filter	1*960 items	VWRA612-6125	3600
10µl Narrow uni Gr tip, 20µl . st, flt	1*960 items	VWRA612-6138	3600
20µl uni Grad tip, Sterile, Filter	1*960 items	VWRA612-6461	3600
50µl uni Grad tip, Sterile, Filter	1*960 items	VWRA612-6467	3600
100µl uni Grad tip, Sterile, Filter	1*960 items	VWRA612-6473	3600
200µl uni Grad tip, Sterile, Filter	1*960 items	VWRA612-6493	3600
200µl Long uni Gr tip, 300ul ., st, flt	1*960 items	VWRA612-6191	3600
1000µl uni Grad tip, Sterile, Filter	1*576 items	VWRA612-6197	2550
1250µl uni Grad tip, Sterile, Filter	1*576 items	VWRA612-6292	2770

STERILE LOW-RETENTION FILTER TIPS

Key features:

- Made of US FDA approved virgin PolyPropylene.
- Especially designed PE filter to prevent cross contamination.
- Racked tips support hands-free tip fitting.
- Each tips are hydrophobic for lowest retention of liquid.

Ordering information:

Sterile Low-retention Filter Bulk

Description	UOM	Cat No.	Price
10µl Long uni Gr tip Sterile, LR, Filter	1*1000 items	VWRA612-6246	2800
10µl Narrow uni Gr tip st, LR, Filter	1*1000 items	VWRA612-6249	2800
20µl uni Grad tip, Sterile, LR, Filter	1*1000 items	VWRA612-6143	2800
50µl uni Grad tip, Sterile, LR, Filter	1*1000 items	VWRA612-6149	2800
100µl uni Grad tip, Sterile, LR, Filter	1*1000 items	VWRA612-6173	2800
200µl uni Grad tip, Sterile, LR, Filter	1*1000 items	VWRA612-6178	2800
200µl Long uni Gr tip, 300ul st, LR, flt	1*1000 items	VWRA612-6262	3293
1000µl uni Grad tip, Sterile, LR, Filter	1*1000 items	VWRA612-6270	3700
1250µl uni Grad tip, Sterile, LR, Filter	1*500 items	VWRA612-6284	1900

Sterile Low-retention Filter Racked

Description	UOM	Cat No.	Price
10µl uni Grad tip, Sterile, LR, Filter	1*960 items	VWRA612-6119	4600
10µl Long uni Gr tip, st, LR, Filter	1*960 items	VWRA612-6247	4600
10µl Narrow uni Gr tip, st, LR, Filter	1*960 items	VWRA612-6260	4600
20µl uni Grad tip, Sterile, LR, Filter	1*960 items	VWRA612-6144	4600
50µl uni Grad tip, Sterile, LR, Filter	1*960 items	VWRA612-6159	4600
100µl uni Grad tip, Sterile, LR, Filter	1*960 items	VWRA612-6174	4600
200µl uni Grad tip, Sterile, LR, Filter	1*960 items	VWRA612-6179	4600
200µl Long uni Gr tip 300ul st LR flt	1*960 items	VWRA612-6263	4600
1000µl uni Grad tip, ST LR, Filter Rack	1*576 items	VWRA612-6272	3021
1250µl uni Grad tip, Sterile, LR, Filter	1*576 items	VWRA612-6293	3500

VWR N95 MASK

POWERED TO PROTECT

Safer


- VWR N95 masks offer proven filtration efficiency with special multiple melt blown layers (>95% for ≥ 0.3-micron droplets)

Softer

- Integrated nose clip for adaptable fitting and comfort
- Latex free, glass fiber free bio compatible material, that is softer yet more protective

Cleaner

- Manufactured in 100K cleanroom environment
- All VWR N95 Masks are pre-sterilized with UV technology



Description	UOM	Cat. No.
N95 Masks	1*25 items	VWRA113-9760

RELOADABLE TIPS:

Key features:

- Readily arrange tips for refilling in the rack
- Time saving & cost saving option
- Available in all variance – filter, non filter, sterile, non sterile and low-retention

Ordering information:

Reloadable Tips

Description	UOM	Cat No.	Price
10µl uni Grad tip	1*960 items	VWRA612-6380	1665
10µl Long uni Grad tip	1*960 items	VWRA612-6406	1665
10µl Narrow uni Grad tip	1*960 items	VWRA612-6434	1665
200µl uni Grad tip	1*960 items	VWRA612-6494	1665
200µl Long uni Grad tip, 300ul .	1*960 items	VWRA612-6522	1829
1000µl uni Grad tip	1*576 items	VWRA612-6544	971
1250µl uni Grad tip	1*576 items	VWRA612-6580	1273

Reloadable Sterile Tips

Description	UOM	Cat No.	Price
10µl uni Grad tip, Sterile	1*960 items	VWRA612-6381	1810
10µl Long uni Grad tip, Sterile	1*960 items	VWRA612-6407	1810
10µl Narrow uni Grad tip, Sterile	1*960 items	VWRA612-6442	1810
200µl uni Grad tip, Sterile	1*960 items	VWRA612-6496	1810
200µl Long uni Gr tip, 300ul ., Sterile	1*960 items	VWRA612-6523	1992
1000µl uni Grad tip, Sterile	1*480 items	VWRA612-6552	1134
1250µl uni Grad tip, Sterile	1*576 items	VWRA612-6222	1470

Reloadable Sterile Low-retention Tips

Description	UOM	Cat No.	Price
10µl uni Grad tip, Sterile, LR	1*960 items	VWRA612-6382	2159
10µl Long uni Grad tip, Sterile, LR	1*960 items	VWRA612-6126	2159
10µl Narrow uni Grad tip, Sterile, LR	1*960 items	VWRA612-6139	2159
200µl uni Grad tip, Sterile, LR	1*960 items	VWRA612-6499	2159
200µl Long uni Gr tip, 300ul ., st, LR	1*960 items	VWRA612-6192	2373
1000µl uni Grad tip, Sterile, LR	1*480 items	VWRA612-6198	1351
1250µl uni Grad tip, Sterile, LR	1*576 items	VWRA612-6294	1726

Individually Wrapped

Description	UOM	Cat No.	Price
200µl uni Grad tip, Sterile (IW)	1*400 items	VWRA612-6506	1795
200µl uni Grad tip, Sterile, LR (IW)	1*400 items	VWRA612-6508	1879
1000µl uni Grad tip, Sterile (IW)	1*200 items	VWRA612-6554	1346
1000µl uni Grad tip, Sterile, LR (IW)	1*200 items	VWRA612-6212	1403

Empty Racks

Description	UOM	Cat No.	Price
10µl tip rack w/o tips	1*10 items	VWRA612-6603	1356
200µl tip rack w/o tips	1*10 items	VWRA612-6604	1356
200XL tip rack w/o tips	1*10 items	VWRA612-6619	1356
1000µl tip rack w/o tips	1*6 items	VWRA612-6627	953
1250µl tip rack w/o tips	1*6 items	VWRA612-6224	1356

Reloadable Filter Tips

Description	UOM	Cat No.	Price
10µl uni Grad tip, Filter	1*960 items	VWRA612-6383	3014
10µl Long uni Grad tip, Filter	1*960 items	VWRA612-6408	3014
10µl Narrow uni Grad tip, Filter	1*960 items	VWRA612-6443	3014
20µl uni Grad tip, Filter	1*960 items	VWRA612-6462	3014
50µl uni Grad tip, Filter	1*960 items	VWRA612-6468	3014
100µl uni Grad tip, Filter	1*960 items	VWRA612-6474	3014
200µl uni Grad tip, Filter	1*960 items	VWRA612-6503	3014
200µl Long uni Grad tip, 300ul ., Filter	1*960 items	VWRA612-6524	3294
1000µl uni Grad tip, Filter	1*480 items	VWRA612-6553	2081
1250µl uni Grad tip, Filter	1*576 items	VWRA612-6223	2297

Reloadable Sterile Filter Tips

Description	UOM	Cat No.	Price
10µl uni Grad tip, Sterile, Filter	1*960 items	VWRA612-6384	3199
10µl Long uni Grad tip, Sterile, Filter	1*960 items	VWRA612-6127	3199
10µl Narrow uni Gr tip, Sterile, Filter	1*960 items	VWRA612-6142	3199
20µl uni Grad tip, Sterile, Filter	1*960 items	VWRA612-6463	3199
50µl uni Grad tip, Sterile, Filter	1*960 items	VWRA612-6469	3199
100µl uni Grad tip, Sterile, Filter	1*960 items	VWRA612-6476	3199
200µl uni Grad tip, Sterile, Filter	1*960 items	VWRA612-6504	3199
200µl Long uni Gr tip, 300ul ., st, flt	1*960 items	VWRA612-6193	3495
1000µl uni Grad tip, Sterile, Filter	1*576 items	VWRA612-6199	2256
1250µl uni Grad tip, Sterile, Filter	1*576 items	VWRA612-6318	2486

Sterile Low Retention Filter Reloadable Tips

Description	UOM	Cat No.	Price
10µl uni Grad tip, Sterile, LR, Filter	1*960 items	VWRA612-6121	3536
10µl Long uni Gr tip, st, LR, Filter	1*960 items	VWRA612-6248	3536
10µl Narrow uni Gr tip, st, LR, Filter	1*960 items	VWRA612-6261	3536
20µl uni Grad tip, Sterile, LR, Filter	1*960 items	VWRA612-6148	3536
50µl uni Grad tip, Sterile, LR, Filter	1*960 items	VWRA612-6172	3536
100µl uni Grad tip, Sterile, LR, Filter	1*960 items	VWRA612-6177	3536
200µl uni Grad tip, Sterile, LR, Filter	1*960 items	VWRA612-6182	3536
200µl Long uni Gr tip 300ul st, LR, flt	1*960 items	VWRA612-6264	3867
1000µl uni Grad tip, Sterile, LR, Filter	1*576 items	VWRA612-6274	2484
1250µl uni Grad tip, Sterile, LR, Filter	1*576 items	VWRA612-6319	2729

Macro Tips

Description	UOM	Cat No.	Price
5000 µl tip for A & VA-FA series pipette	1*100 items	VWRA612-6597	984
10000 µl tip A & VA-FA series pipette	1*100 items	VWRA612-6598	2295
5000µl tip for Smart 5ml pipette series	1*100 items	VWRA612-6599	820
10000 µl tip for Smart 10ml series	1*100 items	VWRA612-6600	1885
10000 µl tip for Smart 10ml series, Sterile	1*100 items	VWRA612-6660	2188

MICRO CENTRIFUGE TUBE

Key features:

0.5 ml, 1.7 ml & 2 ml Graduated Micro Centrifuge Tube

- Manufactured using FDA approved resin
- Certified RNase, DNase, DNA and Pyrogen free
- Free from heavy metals
- Can be centrifuged up to 20,000 RCF
- Withstands temperatures from -80°C to + 121°C
- Easy-seal cap design for smooth single hand operations
- Transparent surface for hindrance-free spectrophotometer reading
- Thin cap membranes in tubes for trouble free needle piercing
- Flat, frosted lid and larger surface on body for easy marking and labelling




Description	UOM	Cat No.	Price
0.5 ml Micro CT, Natural, st	1*1 items	VWRA525-1525-	987
0.5 ml Micro Centrifuge Tube,Amber cl	1*1 items	VWRA525-1527	1074
0.5 ml Micro CT,Amber cl, st	1*1 items	VWRA525-1496	1401
1.7 ml Micro Centrifuge Tube, Natural	1*1 items	VWRA525-1528	987
1.7 ml Micro CT, Natural, st	1*1 items	VWRA525-1529	1197
1.7 ml Micro CT, LR, Sterile	1*1 items	VWRA525-1497	2468
1.7 ml Micro Centrifuge Tube,Amber cl	1*1 items	VWRA525-1531	1179
1.7 ml Micro CT,Amber cl, st	1*1 items	VWRA525-1499	1506
2.0 ml Micro Centrifuge Tube, Natural	1*1 items	VWRA525-1532	599
2.0 ml Micro CT, Natural, st	1*1 items	VWRA525-1533	674
2.0 ml Micro Centrifuge Tube,Amber cl	1*1 items	VWRA525-1535	796
2.0 ml Micro CT,Amber cl, st	1*1 items	VWRA525-1523	963

15 ML & 50 ML GRADUATED CENTRIFUGE TUBE

Key features:

Manufactured using FDA approved polypropylene

- Free from RNase, DNase and Pyrogen
- Free from heavy metals
- Clear transparent surface
- Can be centrifuged up to 17000 RCF & 20000 RCF in 15 ml & 50 ml tube respectively
- Leak proof screw cap
- Tamper proof & chemical resistant graduation printing
- Large writing area for sample identification
- High quality amber colored tubes for highly light sensitive sample ensuring no damage to the samples



Ordering information:

Description	UOM	Cat No.	Price
15 ml Centrifuge Tube with flat Cap, Bulk	1* 500 items	VWRA525-1536	3800
15 ml Centrifuge Tube with flat Cap, Racked	1* 500 items	VWRA525-1537	5000
15 ml CT with flat Cap, Sterile, Bulk Sterile	1* 500 items	VWRA525-1538	4328
15 ml CT with flat Cap, Sterile, Racked	1* 500 items	VWRA525-1539	5600
15 ml CT with flat Cap, Sterile, IW, Individually Wrap	1* 500 items	VWRA525-1540	6646
15 ml CT with flat Cap, Amber cl	1* 500 items	VWRA525-1541	7069
15 ml CT with flat Cap, Amber cl, Racked	1* 500 items	VWRA525-1542	9719
15 ml CT with flat Cap, st, Amber cl	1* 500 items	VWRA525-1543	7731
15 ml CT with flat Cap, st, Amber cl	1* 500 items	VWRA525-1544	9755
15 ml CT with flat Cap, st, IW, Amber cl, Individually Wrapped	1* 500 items	VWRA525-1545	9139
50 ml Centrifuge Tube with flat Cap, Bulk	1* 500 items	VWRA525-1546	5273

Description	UOM	Cat No.	Price
50 ml Centrifuge Tube with flat Cap, Racked	1* 500 items	VWRA525-1547	6885
50 ml CT with flat Cap, Sterile, Bulk	1* 500 items	VWRA525-1548	6590
50 ml CT with flat Cap, Sterile, Racked,	1* 500 items	VWRA525-1549	7869
50 ml CT with flat Cap, Sterile, IW, Individually	1* 500 items	VWRA525-1550	8667
50 ml CT with flat Cap, Amber cl, Bulk	1* 500 items	VWRA525-1551	8561
50 ml CT with flat Cap, Amber cl, Racked	1* 500 items	VWRA525-1552	11103
50 ml CT with flat Cap, st, Amber cl, Bulk	1* 500 items	VWRA525-1553	9341
50 ml CT with flat Cap, st, Amber cl, Racked	1* 500 items	VWRA525-1554	11883
50 ml CT with flat Cap, st, IW, Amber cl,	1* 500 items	VWRA525-1555	11825
50 ml Self Standing CT with flat cap, Bulk	1* 500 items	VWRA525-1556	5432
50 ml Self Standing CT with flat Cap, st, Sterile, Bulk	1* 500 items	VWRA525-1557	6646

TRACE METAL FREE CENTRIFUGE TUBE

Key features:

- Less then 1 PPB (Parts Per Billion) heavy metal content
- Manufactured from ultra pure polypropylene raw material
- No common trace metal ensuring low background in samples
- Traceability of the heavy metal is below ICP detection
- Premium quality PE Cap material (Heavy metal free cap)



Ordering information:

Description	UOM	Cat No.	Price
15 ml CT , Sterile, Heavy metal free cap	1*500 items	VWRA525-1559	5482
15 ml CT , Sterile, IW, Heavy metal free cap	1*500 items	VWRA525-1560	7742
50 ml CT , Sterile, Heavy metal free cap	1*500 items	VWRA525-1562	7936
50 ml CT , Sterile, IW Heavy metal free cap	1*500 items	VWRA525-1563	11148

PCR TUBES

Key features:

Ultra thin walls of 0.35mm thickness ensuring maximum heat transfer for reaction

- Free from PCR inhibitor and Endotoxin
- Certified DNase, RNase & human DNA free
- Specially designed smaller bulk packs to avoid contamination by different users
- Non-reactive internal surface and uniform wall thickness for consistent & error-free results
- Assured evaporation rate of < 5%
- Universal adaptability to thermocyclers of varied brands
- Autoclavable



PCR STRIPS

Key features:

PCR strips come with Snap Caps & 8 tubes strips

- Transparent tube wall for error free results of real time PCR
- Optimized wall thickness for superior heat transfer
- Free from PCR inhibitor, endotoxin, Dnase & RNase,
- Non-reactive internal surface and uniform wall thickness for consistent & error-free results
- Wide range of offering as per need in strips that suits wider applications for PCR and qPCR



Ordering information:

PCR Tubes & PCR Strips:

Description	UOM	Cat No.	Price
0.2ml flat cap clear tubes	1*3000 items	VWRA732-3649	6395
0.2ml PCR 8-Strip Tube Attach Flat Cap	1*120 items	VWRA732-3650	4769
0.2ml PCR 8-Strip Tube Separate Flat Cap	1*120 items	VWRA732-3651	5330
0.2ml PCR 8-Strip Tube with Domed Cap	1*120 items	VWRA732-3652	6115
0.1ml qPCR 8-Strip Tubes Optical Caps	1*120 items	VWRA732-3653	7293
0.1ml qPCR 8-Strip w Tubes Optical Caps	1*120 items	VWRA732-3654	7876
0.2ml qPCR 8-Strip Tubes Optical Caps	1*120 items	VWRA732-3655	7293
0.2ml qPCR 8-Strip w Tubes Optical Caps	1*120 items	VWRA732-3656	7876

PCR PLATES

- Ultra thin walls ensure maximum heat transfer
- Free from PCR inhibitor and Endotoxin
- Certified DNase, RNase & human DNA free
- Specially designed smaller bulk packs to avoid contamination by different users
- Non-reactive internal surface and uniform wall thickness for consistent & error-free results
- Wide range of offering as per need in plates that suits wider applications for PCR and qPCR



TRANSPARENT SEAL

- Strong adhesive prevents evaporation during thermal cycles
- Temperature stability range from -20 to 120° C
- Highly transparent & clear for sample inspection
- Doesn't interfere with the luminescence & fluorescence assay
- Non-pierceable, provides better protection against accidental piercing
- Easy adhesive film doesn't require any pressure to apply on

Ordering information:

Description	UOM	Cat No.	Price
0.2ml, SS Sk, std qPCR 96 Well Pl	1*50 items	VWRA734-3326	12959
0.1ml, SS Sk, Low qPCR 96 Well pl	1*50 items	VWRA734-3327	12959
0.1ml, SS, LP qPCR 96 well pl	1*50 items	VWRA734-3328	12959
PCR 96 well plate sealing Membrane	1*100 items	VWRA734-3330	6171
qPCR 96 well pl optical sealing, Membrane	1*100 items	VWRA734-3331	28387
Aluminum Membrane	1*100 items	VWRA734-3332	6171

RESERVOIR

Key features:

- Made using superior quality polystyrene (PS)
- Convenient pour-off spouts on all four corners.
- Wide base for stability and preventing spillage
- Graduation on inside wall.
- Angled bottom for practically complete withdrawal of sample

Ordering information:

Description	UOM	Cat No.	Price
Reservoir, 10 ML Non sterile Bulk Pack	1*500 items	VWRA613-6725	3900
Reservoir, 10 ML sterile individually wrapped	1*50 items	VWRA613-6726	7158
Reservoir, 10 ML sterile	1*200 items	VWRA613-6727	9391
Reservoir, 25 ML Non Sterile Bulk Pack	1*500 items	VWRA613-6728	4950
Reservoir, 25 ML Sterile individually wrapped	1*50 items	VWRA613-6729	7271
Reservoir, 25 ML Sterile	1*200 items	VWRA613-6730	9391
Reservoir, 50 ML Non Sterile Bulk Pack	1*500 items	VWRA613-6731	6050
Reservoir, 50 ML sterile individually wrapped	1*50 items	VWRA613-6732	7271
Reservoir, 50 ML sterile	1*200 items	VWRA613-6733	9616

LOOP


- Manufactured in clean room on a fully automated production line without human touch
- Designed for dual usage
- Ultra smooth loop surfaces for easy streaking
- Color coded for easy product identification
- Free of lubricants and electrostatic charges for complete liquid transfer
- For spreading liquid samples, smooth rounded surface foot, free of rough edges.
- T-Spreaders are available in two sterile pouch sizes-single wrapped or five spreaders per pack.

Description	UOM	Cat No.	Price
1uL ring and 10uL ring, blue, sterilized, IW	1*500 items	VWRA612-6669	3647
1uL ring and 10uL ring, blue, sterilized	1*1000 items	VWRA612-6670	4488
1uL ring and needle, green, sterilized	1*500 items	VWRA612-6676	3647
1uL ring and needle, green, sterilized	1*1000 items	VWRA612-6677	4488
10uL ring sample and needle, yellow, sterilized, IW	1*500 items	VWRA612-6678	3647
Stick 10uL ring sample and needle, yellow, sterilized	1*1000 items	VWRA612-6679	4488
T-SHAPE wipe for dispersing and wiping, green, sterilized	1*80 items	VWRA612-6680	1122
T-SHAPE wiper for dispersing and wiping, green	1*200 items	VWRA612-6681	2379

VWR BOTTLE TOP DISPENSERS FOR COMPLETE PEACE OF MIND

Key features:

- P: Precise Performance (Manufactured in ISO 8655-5 facility & tested in ISO 17025 accredited lab)
- E: Effortless Dispensing (0.3 grade glass barrel and Teflon 'o'-ring for smooth dispensing experience)
- A: Adaptable Mounting (free adapters to make it fit for almost any bottle)
- C: Contamination free (ease of dismantling and cleaning without disturbing calibration)
- E: Ease of use (360 degree rotating dispensing mouth, double sided graduation)



Ordering information:

Description	UOM	Cat No.	Price
Without Re-Circulation BTD 0.25 to 2.5ml	1*1 items	VWRA612-6628	13340
Without Re-Circulation BTD 0.5 to 5ml	1*1 items	VWRA612-6629	13340
Without Re-Circulation BTD 1 to 10ml	1*1 items	VWRA612-6630	13340
Without Re-Circulation BTD 2.5 to 25ml	1*1 items	VWRA612-6631	19665
Without Re-Circulation BTD 2.5 to 30ml	1*1 items	VWRA612-6632	19665
Without Re-Circulation BTD 5 to 50ml	1*1 items	VWRA612-6633	19665
Without Re-Circulation BTD 5 to 60ml	1*1 items	VWRA612-6634	19665
Without Re-Circulation BTD 10 to 100ml	1*1 items	VWRA612-6667	23460
With Re-Circulation BTD 0.25 to 2.5ml	1*1 items	VWRA612-6635	23690
With Re-Circulation BTD 0.5 to 5ml	1*1 items	VWRA612-6636	23690
With Re-Circulation BTD 1 to 10ml	1*1 items	VWRA612-6637	23690
With Re-Circulation BTD 2.5 to 25ml	1*1 items	VWRA612-6638	33235
With Re-Circulation BTD 2.5 to 30ml	1*1 items	VWRA612-6639	33235
With Re-Circulation BTD 5 to 50ml	1*1 items	VWRA612-6644	33235
With Re-Circulation BTD 5 to 60ml	1*1 items	VWRA612-6649	35190
With Re-Circulation BTD 10 to 100ml	1*1 items	VWRA612-6668	51175