

VWR® SEROLOGY 12 Clinical Centrifuge *BLOOD BANKING CENTRIFUGE*

Part # 77392-164

The VWR® SEROLOGY 12 Clinical Centrifuge is carefully designed to facilitate high-quality blood grouping, typing, cross matching, and other cell-washing procedures. Based on the world's most popular blood banking centrifuge and using the same rotors, the SEROLOGY 12 ensures efficient sample processing, with no individual handling necessary and a convenient automatic pop-lid.

Easy to Use

The fully digital display makes it easy to customize cycle parameters in 5 or 15 second increments. Save up to 10 cycles, enabling easy spinning and simplifying validation. Spin 6 capped or 12 uncapped 75 mm blood tubes with minimal handling (up to 100 mm tubes with additional rotor). Tube positions are numbered for quick identification.

Easy to Monitor

See cycle status at a glance with the VWR® SEROLOGY 12 Clinical Centrifuge's LED lid lighting (patent pending). Bright lights indicate when the centrifuge is ready to load (off), running (on), or done (flashing). Lights continue flashing until the lid is opened.

No Individual Handling

Tubes can be batch filled, centrifuged, transported, incubated, and emptied using the interchangeable heads that also function as holders.

Safe and Durable

The auto-pop lid locks securely during a run, and pops open automatically when the spin is done. The brushless motor and reinforced carbon fiber components ensure years of operation with no routine maintenance.

A Great Value

Unbeatable quality for any budget. The VWR® SEROLOGY 12 Clinical Centrifuge is ready to run with a 2-stage swing-out rotor and 2-year warranty included in the affordable price.



Uses Sero-Fuge™ Rotors

The VWR® SEROLOGY 12 Clinical Centrifuge is compatible with all rotors from the world's most popular blood banking centrifuge.



Sold separately



Included with purchase

Applications	Max Speed	Max Capacity	Included	Warranty
Cell Washing Blood Grouping & Typing Cross Matching	1,040 xg Max RCF 3,600 Max RPM	6 capped or 12 uncapped 75mm tubes*	12 tube rotor <i>*6 place metal rotor (sold separately) required for 100 mm tubes</i>	2 years

VWR® SEROLOGY 12 CLINICAL CENTRIFUGE



Rotor Compatibility

Transition to a new centrifuge without interrupting workflow. The VWR® SEROLOGY 12 Clinical Centrifuge is compatible with all rotors from the world's most popular blood banking centrifuge.

Free 30-Day Evaluation (U.S. Only)

Part Numbers	
Centrifuge	77392-164



6 place metal rotor
0526-001-000 (1 qty)
for 100mm tubes
(sold separately)



12 place rotor
02-001-0-0019 (1 qty)
for 75mm tubes
(1 included)

Specifications	
Sound level	61 dBA
Dimensions (W x D x H)	13 x 15 x 9 in (33 x 38 x 23 cm)
Weight	27 lbs (12.3 kg)
Maximum Speed	3600 RPM
Maximum RCF	1040 xg
Voltage	95-253 VAC
Frequency	50/60 Hz



vwr.com | sales@vwr.com | 1-800-932-5000