

Certificate of Compliance

Product code: ab193971
Product name: Extraction Enhancer Buffer 50X

Lot number: GR3258871-4

Product overview:

Abcam's Extraction Enhancer Buffer 50X is formulated for the preparation of cell or tissues extracts for use in immunoassays. We recommend using ab193971 when having difficulty in pipetting lysate or viscous lysate.

The Extraction Enhancer Buffer 50X should be diluted to 1X in the appropriate buffer before use. If required the Extraction Enhancer Buffer 50X can be used in conjunction of Abcam's Extraction Buffer 5X (ab193970).

The Extraction Enhancer Buffer 50X may precipitate when stored at + 4°C. To dissolve, warm briefly at + 37°C and mix gently.

Species reactivity: Species independent
Storage instructions: Store at +4°C.

Components:

Components	1 ml	5 ml
Extraction Enhancer Buffer 50X	1 x 1ml	1 x 5ml

Abcam warrants that this product conforms to the information contained in this certificate.

The product shelf life is 12 months from date of receipt.

All our products are covered for 12 months by our Abpromise guarantee. For more information, please visit www.abcam.com/abpromise.

Created by (initials): ESTO
Checked by (initials): AGON

Date issued: 01 February 2022

All products are for research use only and are not intended for diagnostic or therapeutic use, not for use in humans.